

rehabhire & sales

YOUR  
LEADING  
HEALTHCARE  
EQUIPMENT  
PARTNER



# VIDA X Wheelchair

MINIMUM  
CHAIR WEIGHT  
**7.5KG**



1300 000 030 | [rehabhire.com.au](http://rehabhire.com.au)

YOUR  
LEADING  
HEALTHCARE  
EQUIPMENT  
PARTNER



The Aspire VIDA X is an ultra-modern, superlightweight transport wheelchair designed for mobility convenience, comfort and style. The wheelchair's ergonomic backrest and seat shape meticulously follow the user's contours, providing maximum comfort and optimal postural support.

With 22 inch sports mag wheels, eye-catching frame and upholstery detail, the VIDA X is a striking modern transport wheelchair solution.

## Lightweight, comfortable, and stylish


- Self-propelled
- Lightweight aluminium body for easy handling and portability
- Enhanced comfort seating via padded seat and backrest with subtle upholstery detail
- Sporty 22 inch mag wheels
- Soft padded armrests
- Swing-away legrests are height adjustable for enhanced flexibility and comfort
- Storage pouch neatly stores essential items
- Choose the matte black frame with striking red upholstery stitching or opt for shiny purple frame with funky purple stitching
- Folds down for easy transport and convenient storage
- Quick release wheels
- Anti-tippers included for added stability
- Lap belt available | WCA848300)

## Specifications

Total Width	690mm
Total Length	1040mm
Total Height	950mm
Seat Width	450mm
Seat Depth	420mm
Anterior Seat Height	500mm
Posterior Seat Height	460mm
Backrest Height	430mm
Armrest Height	200mm
Weight	12.2kg
Transport Weight	7.5kg
SWL	135kg
Warranty	Frame: 5 years   Components: 1 year



**Upholstery**  
Black | Red  
**Frame**  
Matte Black  
MWS449867



**Upholstery**  
Black | Purple  
**Frame**  
Shiny Purple  
MWS449868



Folds for easy transport and storage



Sporty rear wheels



Modern frame colours and upholstery



Storage pouch

Notes

Lined area for notes.

# rehab**hire**&sales

1 300 000 030

320 Lorimer Street Port Melbourne  
contact@rehabhire.com.au  
rehabhire.com.au

# rehab**hire**&sales **REPAIRS AND SERVICE**

REGULAR MAINTENANCE OR REPAIRS

1 300 000 030

Unit 6 | 52 Wirraway Drive Port Melbourne  
service@rehabhire.com.au  
rehabhire.com.au/repairs-and-service

# rehab**installation**

CAN MODIFY EVERY ROOM IN YOUR HOME TO  
ENSURE IT IS ACCESSIBLE, SAFE, AND COMFORTABLE

1 300 000 030

320 Lorimer Street Port Melbourne  
contact@rehabinstallation.com.au  
rehabinstallation.com.au



Registered NDIS Provider