
VICAIR®

POMMEL CUSHION 02

PRESSURE RELIEF AND STABILITY WITH A DETACHABLE MEDIAL THIGH SUPPORT















- Detachable medial thigh support
- Cushion and medial thigh support covered seperately for ease of use

The Vicair Pommel cushion O2 consists of 5 individual compartments including a detachable medial thigh support. This cushion is ideal for individuals with adducted conditions, where the knees or thighs are pulled together. The detachable medial thigh support gently separates the knees and/or thighs relieving the area of pressure that could normally result in pain. The SmartCell filling ensures a softer padded surface between the knees and/or thighs, redistributing pressure in this area and therefore also reducing the risk for pressure injuries. The detachable medial thigh support can be removed to facilitate forward transfers or catheterisation whilst remaining in the wheelchair.

A37A1	LAB	L E C C	VEDC
AVAI	LAB	LE CO	VERS

Comfair Cover and Incotec Cover.

PRODUCT PERFORMANCE			
Skin protection / Anti decubitu	s		
Pressure redistribution	***		
Shear force reduction			
Microclimate control			
Impact & vibration damping	***		
Positioning			
Lateral stability			
Frontal stability	***		
Sliding reduction			
Allows manual adjustments	***		
Increasing endurance			
Maintenance & Safety			
Durability	***		
Reliability			
Functionality when misused			
Easy to maintain			
Level of hygiene			
Functionality & Application			
Comfort	***		
Ultra lightweight	***		
Allows easy transfers			
Available standard sizes	***		
Maximum user weight			

TECHNICAL SPECIFICATIONS			
Thickness	10 cm		
Average weight	850 grams		
Max. user weight	250 kg		
Machine washable	60°C max.		

CUSHION SIZES	
Standard	
On request	

