



SPECIAL MODEL CHART



POMMEL CUSHION O2 - 10 cm

Available sizes:

		WIDTH (WW)											
inch		14	15	16	17	18	19	20	21	22	23	24	
	cm	35	38	40	43	45	48	50	53	55	58	60	
14	35												
15	38												
16	40												
17	43												
18	45												
19	48												
20	50												
21	53												
22	55												
23	58												
24	60												

DEPTH (DD)




Standard Size



Special Size



Custom Size



A = Width

B = Thickness cushion is 10 cm | 4"

C = Depth (Total depth is divided in 1/2 front and 1/2 rear)

D = Thickness pommel is 20 cm | 8"

SIZE = WIDTH (WW) x DEPTH (DD)

rehabhire & sales

320 Lorimer Street Port Melbourne Victoria 3207

t: 03 9646 7200 | e: contact@rehabhire.com.au | rehabhire.com.au