




Bodypoint®

technical bulletin



PART NUMBERS		
SIZE	STANDARD	DYNAMIC
XX-SMALL	SH290XXS	SH295XXS
X-SMALL	SH290XS	SH295XS
SMALL	SH290S	SH295S
MEDIUM	SH290M	SH295M
LARGE	SH290L	SH295L
X-LARGE	SH290XL	SH295XL

TO COMPLETE YOUR PART NUMBER, SELECT THE END-FITTING OPTION BELOW:

-B5



4 - PLASTIC FLAT-MOUNT
END FITTINGS
HKA31-25MM

-B6



4 - PLASTIC CINCH-MOUNT
END FITTINGS
HKA32-25MM

SHEET 1 of 2

note

description

PIVOTFIT HARNESS, MD SERIES, TECHNICAL SPECS

drawing no.

D19-0410-30

date originated

4/ 10/2019

drawn by

D. DOUGAL

rev.

B



Bodypoint®

technical bulletin

FRONT-PULL ORIENTATION

REAR-PULL ORIENTATION

PAD LENGTH


XX-SMALL	24.5cm (9-5/8")
X-SMALL	28cm (11")
SMALL	32.5cm (12-3/4")
MEDIUM	37cm (14-5/8")
LARGE	42.5cm (16-3/4")
X-LARGE	49cm (19-1/4")

LOWER STRAP LENGTH

XX-SMALL	44.5cm (17-1/2")
X-SMALL	
SMALL	54.5cm (21-1/2")
MEDIUM	
LARGE	65cm (25-1/2")
X-LARGE	

COMBINED PAD WIDTH
RANGE

XX-SMALL	11.5 - 19cm (4-1/2 - 7-1/2")
X-SMALL	12.5 - 20.5cm (5 - 8")
SMALL	14 - 24cm (5-1/2 - 9-1/2")
MEDIUM	15 - 25.5cm (6 - 10")
LARGE	17 - 30cm (6-3/4 - 11-3/4")
X-LARGE	19 - 32cm (7-1/2 - 12-1/2")



TOP STRAP ATTACHMENT

SHEET 2 of 2

note

description

PIVOTFIT HARNESS, MD SERIES, TECHNICAL SPECS

drawing no.

D19-0410-30

date originated

4/10/2019

drawn by

D. DOUGAL

rev.

B