



Bodypoint®

technical bulletin





TABLE 1		
SIZE	PART NUMBER	SIZING GUIDE ANKLE CIRCUMFERENCE
X-SMALL	FT240XS	14 - 17cm (5-1/2 - 6-3/4")
SMALL	FT240S	16.5 - 20cm (6-1/2 - 8")
MEDIUM	FT240M	19 - 25.5cm (7-1/2 - 9-1/4")
LARGE	FT240L	23 - 28cm (9 - 11")
X-LARGE	FT240XL	26.5 - 33cm (10-1/2 - 13")



PULL LOOP INSTALLATION / REMOVAL

ADJUSTABLE MOUNTING KIT FT200AS

10-24 x 1/2"
THREAD CUTTING SCREW 8x



FOOTMAN LOOP 4x



PLASTIC SLIDE 4x

note

description

ANKLE HUGGERS TECHNICAL SPECS

drawing no.
D18-0206-30

date originated
2/6/2018

drawn by
D. DOUGAL

rev.
C