



# Slimline Lift Recliner



# Slimline Lift Recliner

Sleek and stylish modern design, the Slimline lift recliner chair puts comfort first and is the perfect addition to any room without being intrusive and bulky. With an individually adjustable headrest, back, lumbar and lift/foot-rest allowing for user personalisation, the chair provides a wide range of comfortable positions. The built-in deep recline feature has a slight tilt in space, allowing users to bring the feet above ninety degrees which is perfect for helping to improve circulation and to relieve pressure. The Slimline comes in three standard sizes with a range of modifications available to give the perfectly fitted chair.



## FEATURES

- 130kg weight capacity
- Adjustable lumbar support
- Controller right side
- Zero gravity recline
- 24 volt battery
- 40mm casters
- Made to measure
- 2 Year Warranty on mechanism, electrics and fabric
- 5 Year Warranty on timber frame
- Available in wide range of Versotela fabrics
- 100% manufactured in Australia by Oscar Furniture

## MODIFICATIONS

CONTROLLER LEFT SIDE

SEAT DEPTH

SEAT WIDTH

ARM HEIGHT

BACK HEIGHT

FOOTREST LENGTH

HEAD AND ARM REST COVERS

SEAT HEIGHT

REMOVABLE SEAT CUSHIONS

CUSHION CUTOUT

MEMORY FOAM

# SPECIFICATIONS



	Seat Height	Seat Foam	Seat Width	Seat Depth	Arm Height	Backrest Height	Footrest Length
MINI STANDARD SIZE							
MIN ADJUSTMENTS				-25mm			-50mm
MEASUREMENT	445mm	80mm	500mm	455mm	165mm	655mm	230mm
MAX ADJUSTMENTS	+30mm	+30mm	+50mm	+50mm	+30mm	+60mm	
A STANDARD SIZE							
MIN ADJUSTMENTS	-30mm	-30mm		-50mm		-30mm	-50mm
MEASUREMENT	475mm	110mm	500mm	495mm	165mm	675mm	230mm
MAX ADJUSTMENTS	+30mm	+30mm	+50mm	+25mm	+30mm	+30mm	
B STANDARD SIZE							
MIN ADJUSTMENTS	-30mm	-30mm		-50mm		-30mm	-50mm
MEASUREMENT	490mm	130mm	500mm	535mm	165mm	695mm	230mm
MAX ADJUSTMENTS			+50mm		+30mm	+30mm	

# FABRIC SELECTION




**VersoTela**  
Fabrics

All VersoTela fabrics come standard with a suite of high-end specifications, such as, antimicrobial coating, mold resistance, scratch resistance and fire resistance.

Our wide range of collections each boast unique properties that are carefully designed to suit every customer's needs. Whether it be built for use in a care facility, or as a functional but stylish addition to the living room, our VersoTela fabric options are sure to tick all the boxes when it comes to crafting your perfect chair.

Please get in touch with one of our friendly sales representatives for fabric samples and to find out more about which fabric is right for you.





Contact us today:  
[sales@oscarfurniture.com.au](mailto:sales@oscarfurniture.com.au)  
03 5381 1404

6 King Drive,  
Horsham VIC 3400

[www.oscarfurniture.com.au](http://www.oscarfurniture.com.au)

