



Order Form: Slimline Lift Chair

Prepared by _____
Client Name _____
Order Date _____
Reference No. _____

Circle required dimensions and configurations

CONFIGURATION

Mechanism	Single Motor Zero-G Mech with adjustable headrest & lumbar		Dual Motor Zero-G Mech with adjustable headrest & lumbar	
Size	MINI		A	B
Wheels	Integrated Retractable Wheels (Standard)		Remove Wheels	
Memory Foam	Seat	Legrest	Backrest	Armrest
Controller Option	Right (Standard)		Left	Double Pocket
Other Options	24V Rechargeable Battery Pack		Scotch Guard	

DIMENSIONS

	Seat Depth (mm)	Seat Height (mm)	Seat Foam (mm)	Seat Width (mm)	Arm Height (mm)	Back Height (mm)	Footrest Length (mm)
Mini Modifications	-25						-50 -25
Mini Standard Size	455	445	80	500	165	655	230
Mini Modifications	+25 +50	+30		+50	+30	+30 +60	
A Modifications	-50 -25	-30				-30	-50 -25
A Standard Size	495	475	110	500	165	675	230
A Modifications	+25	+30		+50	+30	+30	
B Modifications	-50 -25	-30				-30	-50 -25
B Standard Size	535	490	130	500	165	695	230
B Modifications				+50	+30	+30	

NOTE: Reduced seat foam limits Safe working load to 90kg and reduces seat foam warranty to 2 years.

COLOUR AND FINISH

Fabric Inner _____
Outer _____

ADDITIONAL NOTES



www.oscarfurniture.com.au
03 5381 1404
sales@oscarfurniture.com.au

