



B24

Lift Recliner



B24 Lift Recliner

The B24 Lift-Recliner is comfortable for all-day use and can transport the client from room to room without removal from the chair. The 24 volt battery caters for a full day of use. Moving the chair is quick, easy and safe for both the carer and client.



FEATURES

- 130kg weight capacity
- Rear push handle
- 24 volt battery
- Large castors for ease of movement
- Dual motor mechanism
- Made to measure options
- 2 Year Warranty on mechanism, electrics and fabric
- 5 Year Warranty on timber frame
- Available in wide range of Versotela fabrics
- 100% manufactured in Australia by Oscar Furniture

SPECIFICATIONS



	Seat Height	Seat Width	Seat Depth	Arm Height	Backrest Height	Scapular Height	Footrest Length
SIZE A							
MIN ADJUSTMENTS		-150mm	-100mm	-60mm	-120mm	Linked to back height	
MEASUREMENT	476mm	540mm	490mm	180mm	647mm	429mm	N/A
MAX ADJUSTMENTS		+50mm	+75mm	+60mm	+150mm	Linked to back height	+50mm
SIZE B							
MIN ADJUSTMENTS		-150mm	-100mm	-60mm	-150mm	Linked to back height	
MEASUREMENT	496mm	540mm	520mm	190mm	688mm	481mm	N/A
MAX ADJUSTMENTS		+50mm	+50mm	+60mm	+150mm	Linked to back height	+50mm

MODIFICATIONS

- DUAL ACTION
- GLUTERAL RECESS
- LEG REST SIDE WALLS
- ARM HEIGHT
- SEAT DEPTH
- SEAT WIDTH
- ROHO CUTOUT
- BACK HEIGHT
- CONCAVE BACK

FABRIC SELECTION



VersoTela
Fabrics

All VersoTela fabrics come standard with a suite of high-end specifications, such as, antimicrobial coating, mold resistance, scratch resistance and fire resistance.

Our wide range of collections each boast unique properties that are carefully designed to suit every customer's needs. Whether it be built for use in a care facility, or as a functional but stylish addition to the living room, our VersoTela fabric options are sure to tick all the boxes when it comes to crafting your perfect chair.

Please get in touch with one of our friendly sales representatives for fabric samples and to find out more about which fabric is right for you.





Contact us today:
sales@oscarfurniture.com.au
03 5381 1404

6 King Drive,
Horsham VIC 3400

www.oscarfurniture.com.au

