




HEEL PROTECTOR

aspire for... Comfort Safety Prevention

Pressure injuries at the heels are amongst the most commonly acquired pressure injuries in long-term care.

Heels are particularly vulnerable to pressure injury due to direct sustained pressure on the heel from prolonged contact with the supporting surface. As a result, much of the pressure from the weight of the foot and lower leg is concentrated on the small area of skin overlying the heel.

The Aspire Heel Protector is designed to offload pressure from the bony prominence of the heel area, and redistribute it to the larger surface area of the lower leg, reducing the risk of pressure injury at the heel.



1
Year
Warranty



DVA Contracted



Heel Protector is registered in
the ARTG, AUST R ARTG 298451

www.aspirecare.com.au

1300 133 120



aspire HEEL PROTECTOR


- Intelligent leg support zone offloads pressure from the heel and redistributes it to the larger surface area of the lower leg.
- Soft breathable natural cotton fabric outer.
- Cushioning silicone fibre inner promotes airflow for skin microclimate management.
- Seams are located within padding folds to avoid pressure points.
- Adjustable velcro straps offer a snug fit and facilitate one universal size.
- Machine washable.
- Sold as a pair.



PRODUCT DETAILS

Code	Quantity	Sizing
PTA503360	1 x Pair	Universal design fits majority of users

EASY TO APPLY AND REMOVE



ASSOCIATED PRODUCTS

ASPIRE PRESSURE RELIEF CHAIRS



ASPIRE WHEELCHAIR RANGE



AIDACARE ASPIRE LIFTER RANGE



ASPIRE PRESSURE RELIEVING MATTRESSES



www.aspirecare.com.au
 1300 133 120

