

rehabhire & sales

YOUR
LEADING
HEALTHCARE
EQUIPMENT
PARTNER

ELATA[®]

Home Care Mattress

1300 000 030 | rehabhire.com.au

WHY WE EXIST

We believe in zero preventable pressure injuries.

It's ambitious but it's necessary. Pressure injuries cause complications. They also add workload and cost. We act before that pain starts.

WHAT THIS MEANS

- **Care** that preserves dignity.
- **Compassion** that simplifies decisions.
- **Commitment** to better outcomes, every day.

HOW WE DELIVER

Advanced systems that adapt to every patient and pathway. Designed to simplify prescription, enhance comfort, and reduce clinical risk.

We deliver personalised, intuitive pressure care, protecting life with dignity, comfort, and trust.

EVERY PRESSURE INJURY RISK LEVEL. ONE MATCHING SOLUTION.

SIMPLE CHOICES. CLEAR PATHWAYS.

Elata® Range. Quick Category Guide

CATEGORY	MODEL	PURPOSE	PUMP	PG
Home Care	E100, E200	Simple prevention for home and aged care. Default choice in home care pathways.	N/A	9
Clinical Foam	E300	Non-powered. Prevents and treats Stage 1-2 pressure injuries. 9-zone foam with side wall support and heel offloading. Low maintenance.		11
Clinical Hybrid	E400	Use alone for Stage 2 pressure injuries, or step up with a pump to support treatment up to Stage 4. Ideal for home and palliative care. Converts as needs change.	EP10AU	15
Clinical Alternating Hybrid	E500	Alternating air/foam system. Combines comfort with therapy. Treats up to Stage 4 pressure injuries.	EP20AU	19
Clinical Alternating Air	E600	Full alternating air (1-in-2 configuration). Treats up to Stage 4 pressure injuries. Ideal for surgical wards, acute wards, or clients with frailty/reduced mobility.	EP10AU	23
	E700	Advanced alternating (configurable 1-in-4). Auto-setup. Treats up to Stage 4 pressure injuries. For ICU, stable spinal, complex rehab, immobile clients.	EP20AU	27

The EP20AU Pump is fully backwards compatible. It can be used with all Elata pump-ready mattresses (E400 to E700).




Scan to learn more about the range

ELATA[®]

E100

ESSENTIAL FOAM CARE.

Combines comfort foam, high density support foam, and a stable hard base. With side wall reinforcement, cooling fabric, and anti-slip backing, it delivers safe, reliable everyday care for low-risk clients.



S | KS


Cannot be used with a pump
160 kg SWL

ELATA[®]

E200

ENHANCED WITH POCKET SPRINGS.

Builds on the E100 with an integrated pocket spring core beneath its comfort foam. The upgrade delivers cooler airflow, adaptive contouring and greater stability for safer transfers. With the same cooling cover, high density foams, and strong side wall support as the E100, the E200 is the choice for advanced comfort and protection at home.



S | KS

Cannot be used with a pump
160 kg SWL

STANDARD IS JUST THE BEGINNING.

// We are passionate about improving pressure care solutions, delivering advanced foam, hybrid, and alternating air mattresses, along with efficient pumps that set new standards in prevention. Through rigorous research and testing, we aim to exceed industry standards. //
We question every detail and ensure each feature serves a purpose.

Ruben Joyce

Ruben Joyce, CEO and Founder

COMPLETE PEACE OF MIND. LONG COVERAGE. REAL PEOPLE. FAST SUPPORT.

- **10 year conditional pump warranty.** Longest in market.
- **In-house service** with fast diagnostics via USB-C log uploads.
- **Education and tools** that simplify trials and prescriptions.

A pressure care mattress is one intervention in a holistic plan. Always assess the whole patient/client.

Scan to see warranty terms





rehabhire & sales

YOUR
LEADING
HEALTHCARE
EQUIPMENT
PARTNER

1 300 000 030

320 Lorimer Street Port Melbourne
contact@rehabhire.com.au
rehabhire.com.au



rehabhire & sales
REPAIRS AND SERVICE

FOR YOUR
REGULAR
EQUIPMENT
MAINTENANCE
OR REPAIRS

1 300 000 030

Unit 6 | 52 Wirraway Drive Port Melbourne
service@rehabhire.com.au
rehabhire.com.au/repairs-and-service



rehabinstallation

MAKING
EVERYDAY LIFE
IN YOUR HOME
ACCESSIBLE, SAFE,
& COMFORTABLE

1 300 000 030

320 Lorimer Street Port Melbourne
contact@rehabinstallation.com.au
rehabinstallation.com.au