

# aspire REHAB RS *Classic Tilt-In-Space*

The REHAB RS is an ideal solution for the user who requires increased postural support and bodyweight redistribution. A reclined position is easily achieved, making the REHAB RS also a great option for those who fatigue easily.



## INCREASED SUPPORT


The easy-to-adjust, supportive headrest offers the user superior head support, while the gently contoured backrest facilitates posture and positioning.



Infection Control

## EASY TO CLEAN

Wipe clean surfaces ensure infection control compliance.



BACKREST HEIGHT  
ADJUSTABLE

SEAT DEPTH  
ADJUSTABLE

ELEVATING LEGRESTS

## HIGHLY ADJUSTABLE


The wide array of adjustable features fulfil many clinical requirements and are geared towards comfort and support:

- Seat Width
- Seat Depth
- Backrest Height
- Armrest Width & Height
- Elevating Legrests
- Footplate Height
- Wheel & Castor Angle
- Base Frame Length




## ADVANCED COMFORT

Breathable PU coated multi-stretch fabric on seat and back cushions ensure optimal skin microclimate, comfort and improve pressure care outcomes.




26°  
TILT

30°  
RECLINE



## LOW TO MEDIUM LEVEL COMPLEX REHAB SOLUTION

Tilt-In-Space (TIS) functionality promotes increased support and improved positioning. This allows caregivers to position the user for optimal comfort. Adjustable, elevating legrests allow for precise leg positioning.



## SUPPORTIVE & COMFORTABLE

Padded armrests, elevating legrests & a supportive headrest ensure that users of the REHAB RS are supported and comfortable.

- ✓ DVA Contracted
- ✓ AS / NZS Approved
- ✓ Rehab RS is registered in the ARTG, AUST R ARTG 291393



2  
Year  
Warranty  
Frame

1  
Year  
Warranty  
Components

SEAT WIDTH	CODE
410 - 460 mm (16" - 18")	MWS449600
460 - 510 mm (18" - 20")	MWS449610
510 - 560 mm (20" - 22")	MWS449620

Aspire REHAB RS is crash tested to AS/NZS3696.196