

rehabhire & sales

YOUR
LEADING
HEALTHCARE
EQUIPMENT
PARTNER



VIDA Wheelchair

MINIMUM
CHAIR WEIGHT
from
7.5KG



1300 000 030 | rehabhire.com.au

YOUR
LEADING
HEALTHCARE
EQUIPMENT
PARTNER

The Aspire VIDA is an extremely lightweight wheelchair that can be lifted in and out of a vehicle easily by carers. The VIDA range comes in various fun, bold colours and is ideal for quick trips to the shops, appointments, and catching up with friends. The ergonomic backrest and seat shape mimic body contours for comfort and postural support. The high safe working load caters to a wide range of users.

- Attendant or self-propelled
- Three seat width sizes: 400mm, 450mm, and 500mm, each with its own unique colour
- The 400mm seat width is available in mustard
- The 450mm seat width is available in teal and magenta
- The 500mm seat width is available in olive
- Lightweight frame design ensures easy lifting in and out of cars, ideal for carers
- Contoured and padded seating improves comfort and postural support
- Attendant operated hand brakes with parking lock
- Simple and intuitive folding for easy transportation
- Removable leg rests and folding back for compact storage
- Rear storage pocket for essentials
- Quick release wheels
- Anti-tippers included for added stability
- Lap belt available | WCA848300)



Upholstery
Mustard | Black
Frame
Grey



Upholstery
Teal | Black
Frame
Grey



Upholstery
Magenta | Black
Frame
Grey



Upholstery
Olive | Black
Frame
Grey

Attendant Propelled



Self Propelled



Specifications

Attendant Propelled

SKU	MWS449874	MWS449878	MWS449882	MWS449886
Total Width	600mm	650mm	650mm	680mm
Total Length	995mm	995mm	995mm	995mm
Total Height	930mm	930mm	930mm	930mm
Seat Width	400mm	450mm	450mm	450mm
Seat Depth	420mm	420mm	420mm	420mm
Anterior Seat Height	500mm	500mm	500mm	500mm
Posterior Seat Height	460mm	460mm	460mm	460mm
Backrest Height	430mm	430mm	430mm	430mm
Armrest Height	200mm	200mm	200mm	200mm
Rear Wheel Size	13"	13"	13"	13"
Weight	8.8kg	8.8kg	8.8kg	9.3kg
SWL	135kg	135kg	135kg	135kg
Colour	Mustard	Teal	Magenta	Olive
Warranty	Frame: 5 years Components: 1 year			

Self Propelled

SKU	MWS449876	MWS449880	MWS449884	MWS449888
Total Width	630mm	690mm	690mm	710mm
Total Length	1040mm	1040mm	1040mm	1040mm
Total Height	950mm	950mm	950mm	950mm
Seat Width	400mm	450mm	450mm	500mm
Seat Depth	420mm	420mm	420mm	420mm
Anterior Seat Height	500mm	500mm	500mm	500mm
Posterior Seat Height	460mm	460mm	460mm	460mm
Backrest Height	430mm	430mm	430mm	430mm
Armrest Height	200mm	200mm	200mm	200mm
Rear Wheel Size	22"	22"	22"	22"
Weight	7.5kg	7.5kg	7.5kg	8kg
SWL	135kg	135kg	135kg	135kg
Colour	Mustard	Teal	Magenta	Olive
Warranty	Frame: 5 years Components: 1 year			

Notes

rehabhire&sales

1 300 000 030

320 Lorimer Street Port Melbourne
contact@rehabhire.com.au
rehabhire.com.au

rehabhire&sales **REPAIRS AND SERVICE**

REGULAR MAINTENANCE OR REPAIRS

1 300 000 030

Unit 6 | 52 Wirraway Drive Port Melbourne
service@rehabhire.com.au
rehabhire.com.au/repairs-and-service

rehabinstallation

CAN MODIFY EVERY ROOM IN YOUR HOME TO
ENSURE IT IS ACCESSIBLE, SAFE, AND COMFORTABLE

1 300 000 030

320 Lorimer Street Port Melbourne
contact@rehabinstallation.com.au
rehabinstallation.com.au



Registered NDIS Provider