

rehabhire&sales

YOUR
LEADING
HEALTHCARE
EQUIPMENT
PARTNER

AirTec
alternating air system



1300 000 030 | rehabhire.com.au

YOUR
LEADING
HEALTHCARE
EQUIPMENT
PARTNER

The AirTec alternating air system mattresses provide advanced clinical benefits in pressure redistribution and in the treatment of pressure sores, while offering excellent user comfort and support.



AirTec 5 Alternating Pressure Mattress Overlay

The AirTec 5 – 5" overlay provides cost effective care for users with low to medium risk of pressure sores, up to stage II.

Dimensions

Single: 1980mm x 880mm x 130mm
King single: 1980mm x 1050mm x 130mm

Features

- Whisper quiet digital pump
- A & B cells alternating mode - continuous operation
- Static/nursing operation mode
- Alternate or Static seat inflate mode
- Adjustable cycle times
- SWL 160kg

Product Codes

Single: M20-AMP-S
King single: M20-AMP-KS



AirTec 7.5 Alternating Pressure Mattress

The AirTec 7.5 Full mattress replacement features 5" alternating air cells with 2.5" foam underlay. Suitable for users with medium risk/stage III pressure sores.

Dimensions

Single: 1980mm x 880mm x 190mm
King single: 1980mm x 1050mm x 190mm

Features

- Adjustable support weight and comfort settings
- Low pressure alarm with audible buzzer and indicator
- Power failure and service support alarm
- Smart pressure detection technology - prevents the mattress from being over pressurised
- Quick connect plugs
- Full zip-around 2-way stretch nylon TPU cover
- SWL 160kg

Product Codes

Single: M21-AMP-S
King single: M21-AMP-KS

Notice: all products, product specifications and data are subject to change without notice to improve reliability, function, design or otherwise.

Notes

Lined area for notes.

rehab**hire**&sales

1 300 000 030

320 Lorimer Street Port Melbourne
contact@rehabhire.com.au
rehabhire.com.au

rehab**hire**&sales **REPAIRS AND SERVICE**

REGULAR MAINTENANCE OR REPAIRS

1 300 000 030

Unit 6 | 52 Wirraway Drive Port Melbourne
service@rehabhire.com.au
rehabhire.com.au/repairs-and-service

rehab**installation**

CAN MODIFY EVERY ROOM IN YOUR HOME TO
ENSURE IT IS ACCESSIBLE, SAFE, AND COMFORTABLE

1 300 000 030

320 Lorimer Street Port Melbourne
contact@rehabinstallation.com.au
rehabinstallation.com.au



Registered NDIS Provider