

rehabhire&sales

YOUR  
LEADING  
HEALTHCARE  
EQUIPMENT  
PARTNER



# Catalyst 4 & 4C



1300 000 030 | [rehabhire.com.au](http://rehabhire.com.au)

YOUR  
LEADING  
HEALTHCARE  
EQUIPMENT  
PARTNER



**Innovative and Patented Frame Design**

Patented frame concept uses intelligent geometry to create a significantly stronger, lighter and better performing chair.



**Offset Axle Plate**

The Catalyst 4 series vertically adjustable axle plate is mounted forward of the back frame. This small change can significantly improve efficiency.



**Swing Away Footrests**

Enjoy the simplicity and strength of Catalyst's patented swing-away footrest hangers. The tapered engagement keeps the footrests snug and rattle free. Comes with a lifetime warranty.



**Easy to transport. Easy to Push.**

The Catalyst 4 & 4C are the lightest, most durable and easiest to use K4 folders on the market.

**Specifications**

Transport Weight	7.3 kg*
Seat Width	356mm - 559mm
Seat Depth	356mm - 508mm
Front Seat Height	343mm - 533mm
4 Rear Seat Height	318mm - 483mm
4C Rear Seat Height	318mm - 508mm
Top Frame Tubing Size	25.4mm
Bottom Frame Tubing Size	28.6mm
SWL	113.4kg   HD: 158.8kg

\*40cm x 40cm Frame, 127mm x 25.4mm Poly Castors, Height Adjustable T-Arm Receivers. Footrest and Rear Wheels removed.

Notes

Lined area for notes.

# rehabhire&sales

1300 000 030

320 Lorimer Street Port Melbourne  
contact@rehabhire.com.au  
rehabhire.com.au

# rehabhire&sales **REPAIRS AND SERVICE**

REGULAR MAINTENANCE OR REPAIRS

1300 000 030

Unit 6 | 52 Wirraway Drive Port Melbourne  
service@rehabhire.com.au  
rehabhire.com.au/repairs-and-service

# rehabinstallation

CAN MODIFY EVERY ROOM IN YOUR HOME TO  
ENSURE IT IS ACCESSIBLE, SAFE, AND COMFORTABLE

1300 000 030

320 Lorimer Street Port Melbourne  
contact@rehabinstallation.com.au  
rehabinstallation.com.au



Registered NDIS Provider