

rehabhire & sales

YOUR
LEADING
HEALTHCARE
EQUIPMENT
PARTNER



ROGUE²



1300 000 030 | rehabhire.com.au

YOUR
LEADING
HEALTHCARE
EQUIPMENT
PARTNER

Find your limit...and push beyond.

The redefined style and superior design of Rogue2 is the perfect balance of performance, functionality and self-expression. Developed to fit any lifestyle, every element of the chair has been enhanced to provide you with an optimal experience, wherever you go.



Perfectly fit. And perfectly fitting.

The lightest-in-class and engineered to be the perfect fit for your first chair or your next chair. Enhanced strength and durability paired with the right amount of adjustability lets you maximize your efficiency over the life of the wheelchair.





Style and performance.
Redefined.

"It is so important to have a chair that fits me and meets my needs. This is the first wheelchair that feels custom made for me. I can't believe I'm saying this, but I feel so free in it! I'm so happy to be able to go about my day confidently."

- Amanda Steijlen

Founder,
Wheely Big Dreams

Make it your kind of... One of a kind.

Looking good and feeling good go hand in hand. The elegant, sculpted design will make you feel confident wherever you go. Minimalistic or edgy, make your chair one of a kind with the Color Visualiser.

"My passion for outdoor adventure takes me all over the world. Having a chair that can keep up and hold up in remote and rugged locations is essential."

- Joe Stone

Director of Mission,,
Teton Adaptive Sports



YOUR
LEADING
HEALTHCARE
EQUIPMENT
PARTNER

Backrest Design

Adjusting the backrest has never been easier, using one tool and no disassembly, change the backrest angle while the chair is occupied. Relax or adjust your clothing by simply engaging the relaxed position.



Indexed Design

Finally, a camber tube design that makes loading your chair easier. Effortlessly change center of gravity with the precisely aligned Dimple Track. Indexed rear seat height adjustment requires no disassembly.



Wing Design

Optimized wing design lets you get closer for transfers. Corrosion resistant components protect your chair from the elements. Easily get aligned with Ki's intuitive Caster Angle Adjustment System.



Specifications

Transport Weight	4.8kg*
Seat Width	300mm - 510mm
Seat Depth	360mm - 510mm
Front Seat Height	340mm - 520mm
4 Rear Seat Height	340mm - 510mm
Camber	0° 2° 4° 6° 8°
Tubing Size	35mm
SWL	124.7kg

*40cm x 40cm Frame, Short Backrest Height, 101.6mm Poly Castors, Tubular Open Footplate. Rear Wheels removed.



Notes

Lined area for notes.

rehabhire&sales

1 300 000 030

320 Lorimer Street Port Melbourne
contact@rehabhire.com.au
rehabhire.com.au

rehabhire&sales **REPAIRS AND SERVICE**

REGULAR MAINTENANCE OR REPAIRS

1 300 000 030

Unit 6 | 52 Wirraway Drive Port Melbourne
service@rehabhire.com.au
rehabhire.com.au/repairs-and-service

rehabinstallation

CAN MODIFY EVERY ROOM IN YOUR HOME TO
ENSURE IT IS ACCESSIBLE, SAFE, AND COMFORTABLE

1 300 000 030

320 Lorimer Street Port Melbourne
contact@rehabinstallation.com.au
rehabinstallation.com.au



Registered NDIS Provider