

rehabhire&sales

YOUR
LEADING
HEALTHCARE
EQUIPMENT
PARTNER

ottobock.

Avantgarde XXL 2



1300 000 030 | rehabhire.com.au

YOUR
LEADING
HEALTHCARE
EQUIPMENT
PARTNER

The Avantgarde XXL 2 was developed to meet the needs of heavier tall users. Outstanding stability combined with easy operation are its strengths. The robust base structure offers high reliability, and the high quality materials keeps the net weight of the chair as low as possible.

Features and Benefits

- Self and attendant propelled
- Reinforced aluminium tube frame
- Reinforced hollow rims
- Reinforced back support tube with lumbar curve
- Segmented angle adjustable foot support with aluminium footplates



Specifications

SKU	480F45=50000_K
Total Width	695mm - 940mm
Total Length	790mm - 1140mm
Total Height	700mm - 1035mm
Folding Size LxWxH	320mm
Seat Width	460mm – 620mm
Seat Depth	400mm – 560mm
Front Seat Height	440mm – 530mm
Rear Seat Height	410mm – 530mm
Back Height	250mm – 550mm
Lower Leg Length	320mm – 550mm
Front Angle Frame	75°
Approx. Transport Weights	Frame: 11.3kg Drive wheel: >2.15kg Foot Supports: 1.2kg
SWL	180kg
Compliance	ARTG 241133 AS NZS 3695.1:2011
Crash Tested	DIN ISO 7176-19
Warranty	2 years

Notes

Lined area for notes.

rehabhire&sales

1 300 000 030

320 Lorimer Street Port Melbourne
contact@rehabhire.com.au
rehabhire.com.au

rehabhire&sales **REPAIRS AND SERVICE**

REGULAR MAINTENANCE OR REPAIRS

1 300 000 030

Unit 6 | 52 Wirraway Drive Port Melbourne
service@rehabhire.com.au
rehabhire.com.au/repairs-and-service

rehabinstallation

CAN MODIFY EVERY ROOM IN YOUR HOME TO
ENSURE IT IS ACCESSIBLE, SAFE, AND COMFORTABLE

1 300 000 030

320 Lorimer Street Port Melbourne
contact@rehabinstallation.com.au
rehabinstallation.com.au



Registered NDIS Provider