

rehabhire&sales

YOUR
LEADING
HEALTHCARE
EQUIPMENT
PARTNER

A&D
A&D Weighing



Bariatric Stand-on Scale

AMH-400



1300 000 030 | rehabhire.com.au

YOUR
LEADING
HEALTHCARE
EQUIPMENT
PARTNER



General Specifications

- Maximum capacity of 250kg or 400kg
- 100g / 200g resolution provides precise weighing
- BMI function to determine underweight / normal weight / surplus weight
- Easy to read LCD display with white colour LED backlight
- LCD character height: 25mm
- Lock and Hold functions are available as standard
- Platform size is 600mm x 600mm
- Powered by AC adaptor with a rechargeable battery
- Operating time up to 40 hours without backlight, and the charging time is approximately 12 hours.

Technical Specifications

Capacity	400kg	250kg
Resolution	200g	100g

Specifications are subject to change without notice.

Notes

Lined area for notes.

rehabhire&sales

1 300 000 030

320 Lorimer Street Port Melbourne
contact@rehabhire.com.au
rehabhire.com.au

rehabhire&sales **REPAIRS AND SERVICE**

REGULAR MAINTENANCE OR REPAIRS

1 300 000 030

Unit 6 | 52 Wirraway Drive Port Melbourne
service@rehabhire.com.au
rehabhire.com.au/repairs-and-service

rehabinstallation

CAN MODIFY EVERY ROOM IN YOUR HOME TO
ENSURE IT IS ACCESSIBLE, SAFE, AND COMFORTABLE

1 300 000 030

320 Lorimer Street Port Melbourne
contact@rehabinstallation.com.au
rehabinstallation.com.au



Registered NDIS Provider