

rehabhire&sales

YOUR
LEADING
HEALTHCARE
EQUIPMENT
PARTNER

A&D
A&D Weighing



Wheelchair Platform Scale

AMW-300



1300 000 030 | rehabhire.com.au

YOUR
LEADING
HEALTHCARE
EQUIPMENT
PARTNER



General Specifications

- Maximum capacity of 300kg
- 100g resolution provides precise weighing
- Low profile platform with easy wheelchair access from ramps on both sides
- Large steel weighing plate and handrail design, makes it ideal to weigh obese patients
- Easy to read LCD display with white colour LED backlight
- Lock and Hold functions are available as standard
- Handrails and ramps can easily be removed if necessary
- Platform size is 880mm x 840mm
- BMI function to determine underweight/normal weight / surplus weight
- Powered by AC adaptor with a rechargeable battery
- Operating time up to 40 hours without backlight, and the charging time is approximately 12 hours

Technical Specifications

Capacity	300kg
Resolution	100g

Specifications are subject to change without notice.

rehabhire&sales

1 300 000 030

320 Lorimer Street Port Melbourne
contact@rehabhire.com.au
rehabhire.com.au

rehabhire&sales **REPAIRS AND SERVICE**

REGULAR MAINTENANCE OR REPAIRS

1 300 000 030

Unit 6 | 52 Wirraway Drive Port Melbourne
service@rehabhire.com.au
rehabhire.com.au/repairs-and-service

rehabinstallation

CAN MODIFY EVERY ROOM IN YOUR HOME TO
ENSURE IT IS ACCESSIBLE, SAFE, AND COMFORTABLE

1 300 000 030

320 Lorimer Street Port Melbourne
contact@rehabinstallation.com.au
rehabinstallation.com.au



Registered NDIS Provider