

rehabhire&sales

YOUR
LEADING
HEALTHCARE
EQUIPMENT
PARTNER

A&D
A&D Weighing


Bariatric Chair Scale

BCS



1300 000 030 | rehabhire.com.au

YOUR
LEADING
HEALTHCARE
EQUIPMENT
PARTNER



Dimensions

Length	776mm
Width	1027mm
Height	904mm

General Specifications

- High weighing capacity chair scale (up to 500kg)
- Excellent resolution accuracy (100g)
- High quality construction with ample extra large seating
- High quality vinyl padded seating for patient comfort
- Spring loaded foot plate for patient's ease of access
- High strength stainless steel arm rests
- Clear display with standard data output
- Powered with 6 x D cell batteries or optional AC adaptor
- Lockable wheels ensuring patient safety

Technical Specifications

Capacity	500kg
Resolution	100g

Specifications are subject to change without notice.

Notes

Lined area for notes.

rehabhire&sales

1 300 000 030

320 Lorimer Street Port Melbourne
contact@rehabhire.com.au
rehabhire.com.au

rehabhire&sales **REPAIRS AND SERVICE**

REGULAR MAINTENANCE OR REPAIRS

1 300 000 030

Unit 6 | 52 Wirraway Drive Port Melbourne
service@rehabhire.com.au
rehabhire.com.au/repairs-and-service

rehabinstallation

CAN MODIFY EVERY ROOM IN YOUR HOME TO
ENSURE IT IS ACCESSIBLE, SAFE, AND COMFORTABLE

1 300 000 030

320 Lorimer Street Port Melbourne
contact@rehabinstallation.com.au
rehabinstallation.com.au



Registered NDIS Provider