

rehabhire&sales

YOUR
LEADING
HEALTHCARE
EQUIPMENT
PARTNER

A&D
A&D Weighing

Chair Scale

HVL-CS



1300 000 030 | rehabhire.com.au

YOUR
LEADING
HEALTHCARE
EQUIPMENT
PARTNER

Features and Benefits

- Simple operation with 2 button keypad
- Easy to clean chair shell to ensure maximum hygiene
- Long battery life (up to 600hrs) on alkaline cells
- Highly portable with locking castors for patient safety
- High visibility display with easy viewing from multiple angles and auto power off after 3 minutes of inactivity
- Ergonomically designed so minimal effort is required by patient and staff
- Push handles make it convenient for staff when moving the scale
- Australian made



Technical Specifications

Width	750mm
Depth	630mm
Wheels	100mm
Minimum Resolution	20g
Weight Display	Liquid Crystal Display LCD
Chair Shell	White High-strength Moulded Plastic
Base	Grey Die-cast Aluminium Spider
Power	6 x D cell batteries or AC adaptor (optional)
Loadcell	>150% overload
SWL	220kg

Specifications are subject to change without notice.

Notes

Lined area for notes.

rehabhire & sales

1 300 000 030

320 Lorimer Street Port Melbourne
contact@rehabhire.com.au
rehabhire.com.au

rehabhire & sales **REPAIRS AND SERVICE**

REGULAR MAINTENANCE OR REPAIRS

1 300 000 030

Unit 6 | 52 Wirraway Drive Port Melbourne
service@rehabhire.com.au
rehabhire.com.au/repairs-and-service

rehabinstallation

CAN MODIFY EVERY ROOM IN YOUR HOME TO
ENSURE IT IS ACCESSIBLE, SAFE, AND COMFORTABLE

1 300 000 030

320 Lorimer Street Port Melbourne
contact@rehabinstallation.com.au
rehabinstallation.com.au



Registered NDIS Provider