

rehabhire&sales

YOUR
LEADING
HEALTHCARE
EQUIPMENT
PARTNER

A&D
A&D Weighing


Height Measurement Scale

AMP-250



1300 000 030 | rehabhire.com.au

YOUR
LEADING
HEALTHCARE
EQUIPMENT
PARTNER



Features and Benefits

- Maximum capacity of 250kg
- 100g resolution provides precise weighing
- BMI function to determine underweight/normal weight / surplus weight
- Bright LCD display with white colour LED backlight
- Easy to use height measuring stick (880mm-2050mm)
- Non slip base with rubber feet to provide secure footing
- Lock and Hold functions are available as standard.
- Easy to transport with two integrated transport wheels
- Platform size is 370mm x 370mm
- Powered by AC adaptor with a rechargeable battery
- Operating time up to 40 hours without backlight
- Charging time is approximately 12 hours

Technical Specifications

Capacity	250kg
Resolution	100g

Specifications are subject to change without notice.

Notes

Lined area for notes.

rehab**hire**&sales

1 300 000 030

320 Lorimer Street Port Melbourne
contact@rehabhire.com.au
rehabhire.com.au

rehab**hire**&sales **REPAIRS AND SERVICE**

REGULAR MAINTENANCE OR REPAIRS

1 300 000 030

Unit 6 | 52 Wirraway Drive Port Melbourne
service@rehabhire.com.au
rehabhire.com.au/repairs-and-service

rehab**installation**

CAN MODIFY EVERY ROOM IN YOUR HOME TO
ENSURE IT IS ACCESSIBLE, SAFE, AND COMFORTABLE

1 300 000 030

320 Lorimer Street Port Melbourne
contact@rehabinstallation.com.au
rehabinstallation.com.au



Registered NDIS Provider