

rehabhire & sales

YOUR
LEADING
HEALTHCARE
EQUIPMENT
PARTNER

PRIDE[®]
MOBILITY
LIVE YOUR BEST[®]

Jazzy[®]
1450



1300 000 030 | rehabhire.com.au

YOUR
LEADING
HEALTHCARE
EQUIPMENT
PARTNER

Heavy Duty Durability

The Jazzy® 1450 is a heavy-duty, bariatric power chair built for exceptional stability and durability. Engineered to support up to 272 kg, it provides a secure and comfortable ride for users needing a higher weight capacity. With a front-wheel drive system and large, knobby tyres, it delivers smooth performance on even or level terrain, ensuring confident manoeuvrability. Front anti-tip wheels enhance stability, while the rugged frame and efficient motors offer reliable operation for everyday mobility. Designed for strength and comfort, the Jazzy® 1450 is the ideal choice for those seeking a high-capacity power chair for daily use.



Regenerative and Electro-mechanical Braking System



VR2 Joystick Controller



Fold-down Foot Platform

Specifications

Length	970mm	Turning Radius	673mm
Base Width ²	742mm	Battery Requirements ³	2 x 70Ah 12V
Seat Width	609mm	Battery Weight	23.4kg each
Seat Depth	812mm	Range Per Charge Max ¹	<25.2km
Seat to Floor Height	520mm - 570mm	Standard Electronics	90A PG VR2
Drive Train	2 Motor Front Wheel Drive	Specialty Controls	Yes
Drive Wheels	355.6mm Solid	Battery Charger	8A Off-board
Front Castors	101.6mm Solid Anti-tip	Motor Packages	4-pole motors
Rear Castors	228.6mm Solid	Standard Seat Weight 45cm x 45cm	25.8kg
Suspension	Limited	Specialty Seating Options	Synergy® Seating
Maximum Speed ¹	7.2km/h	Weight of Heaviest Piece	88.4kg
Braking System	Regenerative and Electro-mechanical	Available Colours	Red or Blue
Ground Clearance	76mm at Footboard Mount 38mm at Anti-tip Wheel	SWL	273.1kg

1. Varies with user weight, terrain type, battery amp hour (Ah), battery charge, battery condition and tyre condition.

2. Width of base (not including seat).

3. Deep cycle (AGM or gel cell type recommended)

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tyre pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of 3%.

Notes

Lined area for notes.

rehab**hire**&sales

1 300 000 030

320 Lorimer Street Port Melbourne
contact@rehabhire.com.au
rehabhire.com.au

rehab**hire**&sales **REPAIRS AND SERVICE**

REGULAR MAINTENANCE OR REPAIRS

1 300 000 030

Unit 6 | 52 Wirraway Drive Port Melbourne
service@rehabhire.com.au
rehabhire.com.au/repairs-and-service

rehab**installation**

CAN MODIFY EVERY ROOM IN YOUR HOME TO
ENSURE IT IS ACCESSIBLE, SAFE, AND COMFORTABLE

1 300 000 030

320 Lorimer Street Port Melbourne
contact@rehabinstallation.com.au
rehabinstallation.com.au



Registered NDIS Provider