

rehabhire & sales

YOUR
LEADING
HEALTHCARE
EQUIPMENT
PARTNER



Raphael X

Quattro



1300 000 030 | rehabhire.com.au

YOUR
LEADING
HEALTHCARE
EQUIPMENT
PARTNER

Raphael X Quattro Action

4 Motor

Designed for enhanced comfort and customisable support. The Raphael X delivers exceptional comfort with features to suit a diverse range of users. The adjustable backrest includes zipper access for personalised postural adjustment and lateral support, providing improved upright sitting. It's plush design compliments any modern home environment.



Lateral Backrest
Standard



Classic Backrest
Optional

Colour Range



Golden Grain
Microsuede



Marina Grey
Microsuede



Nordic
Stone Vinyl



Custom
Upholstery



Back
Recline



Leg
Raise



Power
Lumbar



Power
Headrest

Functions

Supporting Comfort, Safety, and Independence in Daily Living

The Raphael X range is available in five sizes—Petite, Small, Medium, Large, and Maxi. Each size is designed to provide the right balance of comfort, postural support, and ease of use, ensuring a well-fitted seating experience for a wide range of users.



Adjustable Backrest for Personalised Support

Backrest cushions can be repositioned to accommodate individual postural needs, promoting better alignment and comfort.



Quattro Action Mechanism

Independently adjusts the backrest, leg rest, lumbar support, and headrest to provide tailored comfort for various activities, including reclining, relaxing, and seated support.



Backlit Hand Control with USB Charging

Easy-to-use, illuminated controls improve visibility in low light, while the integrated USB port allows convenient device charging without needing to leave the chair.

- **Built-in Lateral Supports:** Help maintain an upright sitting posture, improving comfort and stability
- **Integrated Footrest Extension:** Helps with fatigue and improves lower limb comfort*
- **Adjustable Seat Depth:** Offers 25mm of depth adjustment, allowing taller users to achieve a better-seated position
- **Protective Covers for High-Wear Areas:** Removable covers for the headrest and armrests to keep your chair looking new
- **Refined, Contemporary Design:** Larger, plush armrests provide extra comfort while maintaining a streamlined look that suits the home environment
- **All-day Battery Backup:** A lithium-ion battery ensures continued operation during power outages, supporting safety and independence
- **Premium Upholstery Options:** Available in modern finishes, including microsuede and leather-look vinyl, with local customisation available by request
- **Classic Backrest Option:** Available as an additional accessory for users who prefer a more traditional support profile

Colour Range

Size	Seat Height	Seat Depth	Seat Width	Back Height	Total Height	Total Depth	Total Width	Distance from Wall	SWL
Petite*	450mm	450mm - 475mm	470mm	630mm	1060mm	840mm	820mm	420mm	130kg
Small	450mm	480mm - 505mm	500mm	630mm	1060mm	870mm	850mm	400mm	130kg
Medium	470mm	500mm - 525mm	525mm	680mm	1100mm	900mm	875mm	400mm	130kg
Large	500mm	530mm - 555mm	550mm	720mm	1100mm	910mm	900mm	420mm	130kg
Maxi	510mm	540mm - 565mm	650mm	675mm	1120mm	910mm	1020mm	420mm	180kg

Warranty: 2 years

NOTE: The petite model is designed without a leg extension for users of shorter stature.

rehabhire&sales

1300 000 030

320 Lorimer Street Port Melbourne
contact@rehabhire.com.au
rehabhire.com.au

rehabhire&sales **REPAIRS AND SERVICE**

REGULAR MAINTENANCE OR REPAIRS

1300 000 030

Unit 6 | 52 Wirraway Drive Port Melbourne
service@rehabhire.com.au
rehabhire.com.au/repairs-and-service

rehabinstallation

CAN MODIFY EVERY ROOM IN YOUR HOME TO
ENSURE IT IS ACCESSIBLE, SAFE, AND COMFORTABLE

1300 000 030

320 Lorimer Street Port Melbourne
contact@rehabinstallation.com.au
rehabinstallation.com.au



Registered NDIS Provider