

rehabhire & sales

2025

YOUR
LEADING
HEALTHCARE
EQUIPMENT
PARTNER



CUSTOM HIRE EQUIPMENT

MANUAL WHEELCHAIRS

Ottobock® Avantgarde 4 DV Package

\$127.00 per week



The Avantgarde lightweight folding wheelchair for active users had been optimised without giving up proven advantages. With numerous individual adjustment options, the Ottobock Avantgarde provides optimum support for your active day to day life.

Size

40cm x 45cm | 45cm x 45cm | 50cm x 45cm

Included in Package

- Breathable tension adjustable backrest
- Desk side panels
- Height-adjustable removable armrests
- Anti-tippers
- Hand rims
- Push handles
- Tip assist

*Rigid back available at additional charge upon request

Crash Tested

AS | NZS 3695.1:2011

SWL

140kg

MANUAL TILT WHEELCHAIRS

Ki Mobility Focus CR Package

\$120.00 per week



The Ki Mobility Focus CR is an innovative tilting wheelchair offering many features to enhance the experience of user and carer.

Size

W: 40cm | 45cm | 50cm

D: 40cm - 50cm

Included in Package

- Ki Mobility Focus CR Base
- Ki Mobility Axiom Posterior Lateral Back
- Ki Mobility 25cm Headrest
- 20" drive wheels with push rims

Crash Tested

AS | NZS 3695.1:2011

Ansi Resna WC-19 Transit Approved

SWL

136kg

POWER WHEELCHAIRS

Quickie® Q300 M Mini Package

\$215.00 per week



Go where others can't with the Q300 M Mini; the narrowest TRUE mid-wheel drive powerchair with a tiny turning circle. Quickie® have taken outdoor performance to the next level with BIG powerchair features to get you further than ever before.

Enjoy a safe and comfortable ride with patented all-wheel suspension that independently cushions every wheel. Climb kerbs of up to 70mm with 30cm drive-wheels. Power through tricky terrain with super-strong 4-pole motors and enjoy an energy-efficient range of up to 30km with 34Ah batteries.

Size

W 45cm | 50cm

D: 40cm - 50cm

Included in Package

- Q300 M Mini Base
- Q300 M Mini Charger
- Q300 M Mini Swing Away Footplates
- Jay J3 Deep Contour Backrest
- Jay Whitmyer Adjust-a-Plush Headrest

Crash Tested

ISO 7176-19/WC19

SWL

136kg



Pride® Mobility Jazzy® 623 2.0 Package

\$250.00 per week



With features including powerful motors that achieve speeds up to 10km/h and improved responsiveness, the Jazzy 623 2.0 is engineered to meet the performance needs of the most active user. The Jazzy 623 2.0 comes standard with TRU-Balance 3, 50° tilt.

Size

40cm | 45cm | 50cm

Included in Package

- Jazzy 623 2.0 Base
- Swing-away or centre mount footplates
- Jazzy 623 2.0 Charger
- ADI Backrest
- ADI Headrest

Crash Tested

AS | NZS 3695.2:2013

SWL

136kg

BARIATRIC POWER WHEELCHAIRS

Pride® Mobility Jazzy® 1450

\$500.00 per week



The Jazzy® 1450 is a heavy-duty, bariatric power chair built for exceptional stability and durability. It provides a secure and comfortable ride for users needing a higher weight capacity. With a front-wheel drive system and large, knobby tyres, it delivers smooth performance on even or level terrain, ensuring confident manoeuvrability. Front anti-tip wheels enhance stability, while the rugged frame and efficient motors offer reliable operation for everyday mobility. Designed for strength and comfort, the Jazzy® 1450 is the ideal choice for those seeking a high-capacity power chair for daily use.

Size

71.12cm

Included in Package

- Pride Mobility Jazzy 1450 base
- Heavy-duty swing away footrests
- Pride Mobility Jazzy 1450 Charger
- Stealth Spectrum Foam SPP Cushion: 71.12cm x 50.8cm
- Stealth Comfort Plus Headrest: 25.4cm
- Padded Backrest
- Lap Belt

SWL

273.1kg

BARIATRIC POWER WHEELCHAIRS

Quantum® Q6 Edge® HD Package

\$460.00 per week



The Q6 Edge HD is engineered to meet the performance needs of the most active user. The Q6 Edge HD accepts a complete range of seating and electronics options to deliver exceptional rehab capability at an excellent value, and is available with our iLevel® technology.

Size

W: 508mm – 609.6mm

D: 50cm - 60cm

Included in Package

- Quantum Q6 Edge HD Base with iLevel
- Swing-away or powered centre mount footplates
- Quantum Q6 Edge HD Charger
- Quantum Tru-Comfort back cushion
- Stealth Spectrum Foam SPP Cushion: 55.88cm x 50.8cm
- Stealth Comfort Plus Headrest
- Attendant controls

* An alternative cushion and or headrest will incur additional hire fees

Crash Tested

AS | NZS 3695.2:2013

SWL

181kg iLevel

Raz-AT Attendant Tilt Shower Commode Chair

\$173.00 per week



The Raz-AT mobile shower commode chair tilts up to 40° to provide extra pressure reduction and comfort. The unique frame design utilizes locking gas springs that help lift the seat to an upright position, making the Raz-AT one of the easiest tilting mobile shower chairs to operate.

The Raz-AT mobile shower commode chair comes standard with the Ischial & Pelvic Alignment System (IPAS), 4" of height adjustment, adjustable-tension back upholstery, molded headrest and fixture, 5" dual-locking casters, flip-up padded armrests, and a commode pan. The easy-to-clean, polished stainless steel frame fits over most toilets and is designed for easy side and front access by the user or attendant.

Size

W: 46cm

D: 40cm - 45cm

Included in Package

- Small head support
- Adjustable back
- Flip-up arm supports
- Lateral thoracic support
- Pelvic belt
- Leg adductors
- BodyPoint Aeromesh calf strap
- Front access moulded seat
- Pan and lid

SWL

160kg



RAZ-AT600 Bariatric Tilt Shower Chair

\$230.00 per week



The Raz-AT600 is a heavy-duty mobile shower commode chair with a stainless steel frame and components. The specially engineered, Weight Assist Spring System (WASS) is available in two weight ranges. It is adjustable to accommodate the user's weight and Center of Mass (CoM) for tilt operation with minimal effort by the attendant.

Size

W: 66cm | 76cm

Included in Package

- Large head support
- Adjustable back
- Flip-up arm supports
- Lateral thoracic support
- Infection control pelvic belt
- Leg adductors
- BodyPoint Aeromesh snap-on calf strap
- Adjustable heel loops
- Front access Visco foam interface
- Pan and lid
- Anti-tippers

SWL

275kg



rehab**hire** & sales

📍 320 Lorimer Street Port Melbourne Victoria 3207

📍 8 Carmart Way Pakenham Victoria 3810

TELEPHONE

1300 000 030

FAX

03 9645 3375

EMAIL

contact@rehabhire.com.au

rehabhire.com.au



Registered NDIS Provider
Authorised supplier for TAC and DVA

TERMS AND CONDITIONS

04/09/2025