

rehabhire&sales

YOUR
LEADING
HEALTHCARE
EQUIPMENT
PARTNER

A&D
A&D Weighing

Blood Pressure Monitor


UA-651SL



1300 000 030 | rehabhire.com.au

YOUR
LEADING
HEALTHCARE
EQUIPMENT
PARTNER

All the essential functions for blood pressure monitoring in a smart compact unit. The UA-651SL was designed with A&D's original concept of blood pressure monitors for easy and accurate monitoring in mind. With a 3-line display and one-touch operation, it retains their fundamental principle of designing easy-to-use monitors. UA-651SL is also equipped with their world leading technologies: IHB and WHO classification indicators. In addition, the latex-free and metal-free SlimFit cuff will provide more comfortable measurement for sensitive skin.



Specifications

Measurement Method	Oscillometric
Measurement Range	Pressure: 0 to 299mmHg Systolic Pressure: 60 to 279mmHg Diastolic Pressure: 40 to 200mmHg Pulse: 40 to 180 beats/minute
Measurement Accuracy	Pressure: ± 3 mmHg Pulse: $\pm 5\%$
Display Type	Digital LCD Pressure/pulse displayed simultaneously
Inflation	Automatic inflation by micro pump Personalised fuzzy logic inflation
Depressurisation	Constant speed exhaust valve
Deflation	Automatic exhaust by solenoid valve
Memory Storage	30
Memory Data Average	Maximum 30 data average of systolic and diastolic values
Operating Conditions	+10°C to +40°C 15 %RH to 85 %RH
Transport Storage Conditions	-20°C to +60°C 10 %RH to 95 %RH
Power Supply	4 x 1.5V batteries (R6P or LR6) or AC adaptor (ECTB-182-O: not included)
Dimensions	Approx. W: 960mm H: 130mm D: 68mm
Unit Weight	Approx. 240g excluding batteries
Cuff	Flat cuff with plastic D-ring Arm circumference: 230mm - 370mm

What is in the Box

- 1 Blood Pressure Monitor
- 1 SlimFit Cuff | 220mm - 320mm
- 1 Carrying Case
- 4 R6P Batteries
- Instruction Manual

Optional Accessories

- Large SlimFit Cuff | 310mm-450mm: UA-CUFBOXLA-C
- AC Adaptor: ECTB-182-O

rehabhire&sales

1 300 000 030

320 Lorimer Street Port Melbourne
contact@rehabhire.com.au
rehabhire.com.au

rehabhire&sales **REPAIRS AND SERVICE**

REGULAR MAINTENANCE OR REPAIRS

1 300 000 030

Unit 6 | 52 Wirraway Drive Port Melbourne
service@rehabhire.com.au
rehabhire.com.au/repairs-and-service

rehabinstallation

CAN MODIFY EVERY ROOM IN YOUR HOME TO
ENSURE IT IS ACCESSIBLE, SAFE, AND COMFORTABLE

1 300 000 030

320 Lorimer Street Port Melbourne
contact@rehabinstallation.com.au
rehabinstallation.com.au



Registered NDIS Provider