

rehabhire&sales

YOUR
LEADING
HEALTHCARE
EQUIPMENT
PARTNER

A&D
A&D Weighing

Professional Scale

PW-200-FG



1300 000 030 | rehabhire.com.au


YOUR
LEADING
HEALTHCARE
EQUIPMENT
PARTNER



General Specifications

- High capacity platform scale (up to 200kg)
- Excellent resolution accuracy (20g)
- Large non-slip weigh platform
- Powered by AC adaptor or 4 x C cell batteries
- Easy Portability
- Clear LCD display, and easy to use keyboard
- RS-232C interface option available
- Made in Australia

Dimensions | mm



Technical Specifications

Capacity	200kg
Resolution	20g

Specifications are subject to change without notice.

rehabhire & sales

1 300 000 030

320 Lorimer Street Port Melbourne
contact@rehabhire.com.au
rehabhire.com.au

rehabhire & sales **REPAIRS AND SERVICE**

REGULAR MAINTENANCE OR REPAIRS

1 300 000 030

Unit 6 | 52 Wirraway Drive Port Melbourne
service@rehabhire.com.au
rehabhire.com.au/repairs-and-service

rehabinstallation

CAN MODIFY EVERY ROOM IN YOUR HOME TO
ENSURE IT IS ACCESSIBLE, SAFE, AND COMFORTABLE

1 300 000 030

320 Lorimer Street Port Melbourne
contact@rehabinstallation.com.au
rehabinstallation.com.au



Registered NDIS Provider