

rehabhire & sales



YOUR  
LEADING  
HEALTHCARE  
EQUIPMENT  
PARTNER

# OSLO Seat Walker



1300 000 030 | [rehabhire.com.au](http://rehabhire.com.au)

YOUR  
LEADING  
HEALTHCARE  
EQUIPMENT  
PARTNER



Sleek lines, elegant curves and eye-catching shades give the OSLO seat walker sophisticated poise. The aero-space grade aluminium frame is highly durable and ultra lightweight for easy manoeuvrability, folding and lifting, while tool-free adjustments can suit individual user needs.

- Constructed with aerospace-grade aluminium for exceptional strength, the 6.25kg minimum weight makes it easy to lift and transport
- Tool-free handle and backrest height adjustment accommodates users of varying stature, with push-button locking for secure customisation
- Contoured, non-slip hand grips to reduce wrist strain and adjustable backrest for tailored support
- Puncture-resistant PU wheels with swivel front castors and fixed rear wheels improve manoeuvrability
- One-click folding mechanism collapses the rollator into a slim profile (20cm width) for effortless car boot storage or travel
- Secure zippered storage bag with 5kg capacity and built-in cane holder






**Fold**  
Pull handle up to collapse & lock



**Adjust**  
Press and hold to adjust



**Brake**  
Lock: down  
Unlock: up

## Specifications

	Total Length	Total Width	Total Height	Seat Width	Seat Depth	Seat Height	Front Wheel	Rear Wheel	Folded Width	Frame Weight	Total Weight	SWL
● WAF705320PU	660mm	540mm	830mm - 940mm	450mm	205mm	585mm	230mm	195mm	200mm	6.25kg	6.9kg	136kg
● WAF705320GR												

Warranty

Frame: 2 years | Components: 12 months

Notes

Lined area for notes.

# rehabhire&sales

1300 000 030

320 Lorimer Street Port Melbourne  
contact@rehabhire.com.au  
rehabhire.com.au

# rehabhire&sales **REPAIRS AND SERVICE**

REGULAR MAINTENANCE OR REPAIRS

1300 000 030

Unit 6 | 52 Wirraway Drive Port Melbourne  
service@rehabhire.com.au  
rehabhire.com.au/repairs-and-service

# rehabinstallation

CAN MODIFY EVERY ROOM IN YOUR HOME TO  
ENSURE IT IS ACCESSIBLE, SAFE, AND COMFORTABLE

1300 000 030

320 Lorimer Street Port Melbourne  
contact@rehabinstallation.com.au  
rehabinstallation.com.au



Registered NDIS Provider