

rehabhire & sales

YOUR
LEADING
HEALTHCARE
EQUIPMENT
PARTNER



Wheelchair Accessories



1300 000 030 | rehabhire.com.au

YOUR
LEADING
HEALTHCARE
EQUIPMENT
PARTNER



Anti-tip Bars

Ensures chair does not tilt over when traveling across inclines and uneven surfaces.

AVAILABLE FOR

Vida • Vida X • Vida MG • Socialite • Evoke 2 • Evoke 2 JNR • Transit 3 • Assist 3 • Assist 3 Ward • Rehab RX2 • Rehab RX2 JNR



Attendant Brake Kits

Transit conversion Kit includes 12" wheels and attendant brakes , or attendant brakes kit without wheels.

AVAILABLE FOR

Evoke 2



Safety Lap Belt

Suitable for supporting users with low tone or weakness. Metal push button for secure fit.

AVAILABLE FOR

Vida • Vida X • Vida MG • Socialite • Evoke 2 • Evoke 2 JNR • Transit 3 • Assist 3 • Assist 3 Ward • Rehab RX2 • Rehab RX2 JNR



Contoured Padded Overlay

Increased comfort, ideal for users in their wheelchair for longer periods of time. Breathable washable cover.

AVAILABLE FOR

Evoke 2 • Evoke 2 JNR • Transit 3 • Assist 3 • Assist 3 Ward



Drum Brake Kit

This upgrade kit can be retrofitted, providing the attendant with more control on sloped surfaces.

AVAILABLE FOR

Evoke 2 • Evoke 2 JNR • Rehab RX2 • Rehab RX2 JNR



Oxygen Bottle Holder

Holds oxygen bottle in place. Designed to hold C size bottles.

AVAILABLE FOR

Evoke 2 • Evoke 2 JNR • Rehab RX2 • Rehab RX2 JNR • Transit 3 • Assist 3



Stump Support

For amputee support.

AVAILABLE FOR

Evoke 2 • Evoke 2 JNR • Rehab RX2 • Rehab RX2 JNR • Transit 3 • Assist 3



Elevating Legrest

Padded legrests for improved comfort.

AVAILABLE FOR

Evoke 2 • Evoke 2 JNR • Transit 3 • Assist 3 • Assist 3 Ward • Rehab RX2 • Rehab RX2 JNR



IV Pole

Facilitates intravenous infusion.

AVAILABLE FOR

Evoke 2 • Evoke 2 JNR • Transit 3 • Assist 3 • Assist 3 Ward • Rehab RX2 • Rehab RX2 JNR



Extension Lever

Extends brake lever for user accessibility.

AVAILABLE FOR

Evoke 2 • Evoke 2 JNR • Transit 3 • Assist 3 • Assist 3 Ward



Calf Strap

Provides extra support for user's calves.

AVAILABLE FOR

Evoke 2 • Evoke 2 JNR • Transit 3 • Assist 3 • Assist 3 Ward • Rehab RX2 • Rehab RX2 JNR



Base Frame Extension Kit

Extends stump support to cater for taller users.

AVAILABLE FOR

Rehab RX2 • Rehab RX2 JNR



Depth & Height Adjustable Armrests

Provides multiple support positions for users arms.

AVAILABLE FOR

Evoke 2



Extended Back Canes

Extends back canes up to 50mm to facilitate taller carers.

AVAILABLE FOR

Evoke 2 • Evoke 2 JNR



Stroller Handle

Angle adjustable handle allows greater ease of use for a carer.

AVAILABLE FOR

Rehab RX2 • Rehab RX2 JNR



Lateral Chest Support

Facilitates postural positioning.

AVAILABLE FOR

Evoke 2 • Evoke 2 JNR • Rehab RX2



Headrest Assembly

Deluxe angle adjustable head support. Helps to promote user posture and comfort.

AVAILABLE FOR

Evoke 2



Gutter Armrest Pad

Provides extra support for users with weak or unsupported arms.

AVAILABLE FOR

Evoke 2



Adjustment Kit

Allows for easy, tool-free adjustment of a wheelchair.

AVAILABLE FOR

Rehab RX2 • Rehab RX2 JNR



Meal Tray

Provides sturdy surface for meal time.

AVAILABLE FOR

Rehab RX2



Mobility Accessory Bags

Keep personal items close and secure when you're on the move.

AVAILABLE FOR

Socialite • Vida • Vida X • Vida MG • Evoke 2 • Lite • Transit 3 • Assist 3



Leopard - WAA696500



Plaid - WAA696520



Floral - WAA696510



Blue - WAA696530

Aspire Rehab RX2 Wheelchair Accessories

○ Included as Standard ● Compatible

Code	Description	Rehab RX2	Rehab RX2 JNR
Stump Supports			
WCA848220	Left - 300/320/350/400/410mm	●	●
WCA848230	Right - 300/320/350/400/410mm	●	●
WCA848252	Left/Right - 410 to 460mm	●	●
WCA848254	Left/Right - 510 to 610mm	●	●
Lateral Chest Support Assembly			
WCA848450	Left	●	
WCA848460	Right	●	
Elevating Legrests			
WCA848100	Left - 300/320/350/400/410mm	●	●
WCA848110	Right - 300/320/350/400/410mm	●	●
WCA848120	Left - 450/460/500/510mm	●	
WCA848130	Right - 450/460/500/510mm	●	
Oxygen Bottle Holders			
WAA695600	Oxygen Bottle Bag	●	●

Code	Description	Rehab RX2	Rehab RX2 JNR
Rear Wheels with Drum Brake Kits			
WCA848510	16"	●	●
WCA848520	20"	●	○
WCA848500	24"	○	
Base Frame Extension Kits			
WCA848530	50mm - Black	●	
WCA848540	100mm - Black	●	
Other			
WCA848550	Tool-free Adjustment Kit	●	●
WCA848570	Meal Tray - 460mm to 510mm	●	
WCA848270	Anti-tip Bars	○	○
WCA848280	IV Pole	●	●
WCA848300	Safety Lap Belt - Metal Seatbelt Buckle	●	●
WCA848580	Calf Strap	●	●
WCA848590	Stroller Handle	●	●

*Note Side transfers not possible when gutter armrests are attached

Aspire Vida Range Wheelchair Accessories

Code	Description
WCA848275	Anti-tip Bars - Pair
WCA848300	Safety Lap Belt - Metal Seatbelt Buckle

Aspire Socialite Wheelchair Accessories

Code	Description
WCA848270	Anti-tip Bars - Pair

Aspire Evoke 2, Transit 3 & Assist 3 Wheelchair Accessories

Code	Description	Evoke 2	Evoke 2 JNR	Transit 3	Assist 3
Contoured Padded Overlays					
WCA848325	300mm		●		
WCA848330	350mm	●			
WCA848340	400mm	●		●	●
WCA848350	450mm	●		●	●
WCA848360	500mm	●		●	●
WCA848370	570mm	●		●	●
WCA848380	610mm	●		●	●
Stump Supports					
WCA848220	Left - 300 to 410mm	●	●		
WCA848230	Right - 300 to 410mm	●	●		
WCA848252	Left/Right - 410 to 510mm	●			
WCA848254	Left/Right - 510 to 610mm	●			
WCA848260	Extension Bracket	●	●		
Elevating Legrests					
WCA848100	Left - 300 to 410mm	●	●		
WCA500080	Left - 400mm			●	●
WCA848110	Right - 300 to 410mm	●	●		
WCA500085	Right - 400mm			●	●
WCA848120	Left - 450 to 510mm	●			
WCA500070	Left - 450 to 610mm			●	●
WCA848130	Right - 450 to 510mm	●			
WCA500075	Right - 450 to 610mm			●	●
WCA848140	Left - 570/610mm	●			
WCA848150	Right - 570/610mm	●			
WCA848080	Left - Evoke 2 JNR		●		
WCA848090	Right - Evoke 2 JNR		●		
Lateral Chest Support Assembly					
WCA848430	Left	●	●		
WCA848440	Right	●	●		
Depth & Height Adjustable Armrests					
WCA848390	Left	●			
WCA848400	Right	●			

Code	Description	Evoke 2	Evoke 2 JNR	Transit 3	Assist 3
Gutter Armrest Pads*					
WCA848410	Left	●			
WCA848420	Right	●			
Rear Wheels with Drum Brake Kits					
WCA848510	16"		●		
WCA848520	20"		●		
WCA848500	24"	●			
Brake Kits					
WCA848592	Attendant Brakes Kit (Brakes Only)	●			
	Transit Conversion Kit	●			
WCA848594	(Brakes & 12" Wheels)				
Brake Extension Levers					
WCA848310	Left	●	●	●	●
WCA848320	Right	●	●	●	●
Anti-tip Bar					
WCA848270	Anti-tip Bars - Pair	●	●		
WCA500050	Anti-tip Bars - Pair			●	●
WCA848272	Depth -adjustable Anti-tip Bars - Pair	●			
Other					
WCA848290	Oxygen Bottle Holder	●	●		
WCA848280	IV Pole	●	●		
WCA848285	IV Pole with Mounting Clamp			●	●
WCA848300	Safety Lap Belt - Metal Seatbelt Buckle	●	●		
WCA500055	Safety Lap Belt - Metal Seatbelt Buckle			●	●
WCA848470	Extended Back Cane - Pair - 19"/21"	●	●		
WCA848480	Headrest Assembly 400 to 450mm	●	●		
WCA848481	Headrest Assembly 500 to 610mm	●			
WCA848580	Calf Strap	●	●		
WCA500065	Calf Strap			●	●
WCA500060	Tool Kit			●	●
AAL848590	Security Mod			●	●

*Side transfers not possible when gutter armrests are attached



Notes

Lined area for notes.

rehabhire&sales

1 300 000 030

320 Lorimer Street Port Melbourne
contact@rehabhire.com.au
rehabhire.com.au

rehabhire&sales **REPAIRS AND SERVICE**

REGULAR MAINTENANCE OR REPAIRS

1 300 000 030

Unit 6 | 52 Wirraway Drive Port Melbourne
service@rehabhire.com.au
rehabhire.com.au/repairs-and-service

rehabinstallation

CAN MODIFY EVERY ROOM IN YOUR HOME TO
ENSURE IT IS ACCESSIBLE, SAFE, AND COMFORTABLE

1 300 000 030

320 Lorimer Street Port Melbourne
contact@rehabinstallation.com.au
rehabinstallation.com.au



Registered NDIS Provider