

rehabhire & sales

YOUR  
LEADING  
HEALTHCARE  
EQUIPMENT  
PARTNER



# Ninja Carbon Fibre Seat Walker



1300 000 030 | [rehabhire.com.au](http://rehabhire.com.au)

YOUR  
LEADING  
HEALTHCARE  
EQUIPMENT  
PARTNER




The Aspire Ninja is a lightweight seat walker designed for everyday use. Made from carbon fibre, it's easy to lift in and out of a car and has large front wheels for greater mobility outdoors.



- Rear wheels are positioned for maximum clearance during entry and exit from the frame
- Folds easily for convenient transport and compact storage
- The large 254mm wheels provide stability and smooth movement across a variety of terrains
- Soft-Touch Brakes
- Lock the rear wheels with soft touch brakes to enhance safety and stability during use
- Ergonomic Handles
- The ergonomic flat top handles provide a comfortable grip during extended use and a stable resting place for arms while seated
- Adjustable handles allow for personalised positioning, helping to maintain proper posture
- Utility bag holds up to 5kg, offering convenient storage for personal items
- Compact design improves turning, offering smoother navigation in tight spaces




## Accessories




Adjustable  
Height IV Pole  
WAA695650



Oxygen Bottle  
Holder Bag  
WAA695600



Night Light LED  
Lamp  
WAA695630



Crutch & Cane  
Holder  
Included

## Specifications

	Total Length	Total Width	Total Height	Seat Width	Seat Height	Front Wheel	Rear Wheel	Folded Width	Frame Weight	Total Weight	SWL
Medium: WAF705380	710mm	625mm	860mm - 945mm	460mm	540mm	254mm	203mm	245mm	5.4kg	6.35kg	150kg
Tall: WAF705390	710mm	625mm	860mm - 1025mm	460mm	620mm	254mm	203mm	245mm	5.5kg	6.4kg	150kg

Compliance

ARTG 271255

Warranty

Frame: 5 years | Components: 12 months

Notes

Lined area for notes.

# rehabhire&sales

1300 000 030

320 Lorimer Street Port Melbourne  
contact@rehabhire.com.au  
rehabhire.com.au

# rehabhire&sales **REPAIRS AND SERVICE**

REGULAR MAINTENANCE OR REPAIRS

1300 000 030

Unit 6 | 52 Wirraway Drive Port Melbourne  
service@rehabhire.com.au  
rehabhire.com.au/repairs-and-service

# rehabinstallation

CAN MODIFY EVERY ROOM IN YOUR HOME TO  
ENSURE IT IS ACCESSIBLE, SAFE, AND COMFORTABLE

1300 000 030

320 Lorimer Street Port Melbourne  
contact@rehabinstallation.com.au  
rehabinstallation.com.au



Registered NDIS Provider