

rehabhire & sales

YOUR
LEADING
HEALTHCARE
EQUIPMENT
PARTNER



Rehab RX2

Tilt-in-Space



1300 000 030 | rehabhire.com.au

YOUR
LEADING
HEALTHCARE
EQUIPMENT
PARTNER


The Rehab RX2 Tilt-in-Space Wheelchair is designed to support individuals with complex seating, positioning, or pressure care needs. Its highly adjustable features and compatibility with pressure care cushions enable clinicians to tailor the chair to each user's unique requirements, facilitating effective posture management, enhanced comfort, and more efficient daily care routines.




Features and Benefits

- **Tilt-in-space up to 35°:** Allows postural variation throughout the day, facilitating pressure redistribution and enhancing comfort during prolonged seated periods
- **Lightweight aluminium frame** (from 27.7kg): Minimises carer effort during pushing or manual handling, while supporting safer daily mobility and transport requirements
- **Customisable seating dimensions:** Adjustable seat width, depth, backrest height, armrest position, and headrest placement allow for tailored postural support to meet individual user needs
- **Contoured cushions with multi-stretch covers:** Provide stable support and help reduce pressure build-up, promoting skin protection and overall user comfort
- **Elevating and swing-away legrests:** Assist with leg positioning to improve blood circulation, support foot and ankle posture, and enable easier side or front transfers
- **Removable rear wheels and height-adjustable castors:** Enhance manoeuvrability, simplify storage in tight spaces, and support vehicle loading
- **Brake options for user and carer:** Dual systems, including push-to-lock and drum brakes, provide users with independence and carers with added safety control
- **Transport ready:** Crash-tested to AS/NZS 3691.19 standards for use in vehicle transport with appropriate restraint systems
- **Anti-tippers included as standard:** Increase rear stability when tilting or navigating slopes and ramps


Adjustments




Backrest Height



Seat & Backrest Depth



Elevating Legrest



Tilt-in-space
35° Tilt | 30° Recline

Seat Width	Seat Depth	Seat Height	Total Width	Backrest Height	Armrest Height	Weight	SWL
410mm - 460mm	450mm - 580mm	430mm - 505mm	570mm - 640mm	410mm - 730mm	170mm - 300mm	27.1kg	150kg
460mm - 510mm	450mm - 580mm	430mm - 505mm	620mm - 690mm	410mm - 730mm	170mm - 300mm	30kg	150kg
510mm - 560mm	450mm - 580mm	430mm - 505mm	670mm - 740mm	410mm - 730mm	170mm - 300mm	32.3kg	150kg
Warranty			Frame: 2 years Components: 12 months				

Notes

Lined area for notes.

rehab**hire**&sales

1 300 000 030

320 Lorimer Street Port Melbourne
contact@rehabhire.com.au
rehabhire.com.au

rehab**hire**&sales **REPAIRS AND SERVICE**

REGULAR MAINTENANCE OR REPAIRS

1 300 000 030

Unit 6 | 52 Wirraway Drive Port Melbourne
service@rehabhire.com.au
rehabhire.com.au/repairs-and-service

rehab**installation**

CAN MODIFY EVERY ROOM IN YOUR HOME TO
ENSURE IT IS ACCESSIBLE, SAFE, AND COMFORTABLE

1 300 000 030

320 Lorimer Street Port Melbourne
contact@rehabinstallation.com.au
rehabinstallation.com.au



Registered NDIS Provider