

rehabhire & sales

YOUR
LEADING
HEALTHCARE
EQUIPMENT
PARTNER



Assist-A-Tray



1300 000 030 | rehabhire.com.au

YOUR
LEADING
HEALTHCARE
EQUIPMENT
PARTNER

Features and Benefits

- **Ergonomic Safety Handle:** Makes standing and sitting easy and safe.
- **Swivel Tray:** Multi-use laptop/TV tray with cup holder and utensil compartment. Tray swivels 360°.
- **Adaptable:** Adjusts in length and height to fit either side of your favorite couch or chair.
- **Super Grip Rubber Pads:** Ensure stability and protect floor from damage.
- **Simple Assembly:** No additional tools required



Easy Stand Safety Handle



Adjustable Height and Base Depth



Arm Pivots and Tray Swivels

Specifications

| | |
|---------------|--------------------------------------|
| Tray Height | 660.4mm - 812.8mm |
| Handle Height | 863.6mm - 1016mm |
| Tray Size | 406.4mm x 457.2mm |
| Handle Size | H: 152.4mm W: 152.4mm D: 31.75mm |
| Base Length | 508mm - 914.4mm |
| Base Pads | 101.6mm - 177.8mm |
| Weight | 11.343kg |
| SWL | Handle: 113.4kg Tray: 13.61kg |
| Warranty | |

In the Box



rehabhire&sales

1300 000 030

320 Lorimer Street Port Melbourne
contact@rehabhire.com.au
rehabhire.com.au

rehabhire&sales **REPAIRS AND SERVICE**

REGULAR MAINTENANCE OR REPAIRS

1300 000 030

Unit 6 | 52 Wirraway Drive Port Melbourne
service@rehabhire.com.au
rehabhire.com.au/repairs-and-service

rehabinstallation

CAN MODIFY EVERY ROOM IN YOUR HOME TO
ENSURE IT IS ACCESSIBLE, SAFE, AND COMFORTABLE

1300 000 030

320 Lorimer Street Port Melbourne
contact@rehabinstallation.com.au
rehabinstallation.com.au



Registered NDIS Provider