

rehabhire & sales

YOUR
LEADING
HEALTHCARE
EQUIPMENT
PARTNER

TopGun
Mobility
Mobility & Homecare Specialists

Athena Lift Recliner



1300 000 030 | rehabhire.com.au

YOUR
LEADING
HEALTHCARE
EQUIPMENT
PARTNER

Indulge in Supreme Comfort!

Experience ultimate relaxation and comfort with the luxurious Athena Recliner Lounge. Sink into plush cushioning as you recline, supported by ergonomically designed features that cradle your body in blissful comfort.

Features and Benefits

- OKIN motor
- 35D high density memory foam with spring packet, S-spring on seat and backrest
- USB charged remote with backlight
- Suede look water, oil, and dust resistant tech fabric
- Washable protective headrest/armrest covers included
- Side pocket
- Battery Backup

Specifications

Total Height	1030mm
Total Width	860mm
Total Depth	930mm
Seat Width	480mm
Seat Depth	510mm
Colours	Grey Light Grey
Weight	45kg
SWL	150kg

Note:

- All specifications subject to change without notice.
- The Information contained herein is correct at the time of publication; due to our commitment constantly improve products during development, we reserve right alter specifications without prior notice.
- The product images shown are for illustration purposes only and may not be an exact representation of the product.
- It is the consumers responsibility to ensure they confirm the accuracy of specifications listed here to the product they are purchasing.
- Speed and range vary with user weight, terrain type, battery charge, battery condition and tyre pressure.



Grey



Light Grey



Notes

Lined area for notes.

rehabhire&sales

1300 000 030

320 Lorimer Street Port Melbourne
contact@rehabhire.com.au
rehabhire.com.au

rehabhire&sales **REPAIRS AND SERVICE**

REGULAR MAINTENANCE OR REPAIRS

1300 000 030

Unit 6 | 52 Wirraway Drive Port Melbourne
service@rehabhire.com.au
rehabhire.com.au/repairs-and-service

rehabinstallation

CAN MODIFY EVERY ROOM IN YOUR HOME TO
ENSURE IT IS ACCESSIBLE, SAFE, AND COMFORTABLE

1300 000 030

320 Lorimer Street Port Melbourne
contact@rehabinstallation.com.au
rehabinstallation.com.au



Registered NDIS Provider