

rehabhire & sales

YOUR
LEADING
HEALTHCARE
EQUIPMENT
PARTNER

BetterLiving
by novis®

The Carbon

Strength, style, and ultralight design



1300 000 030 | rehabhire.com.au

YOUR
LEADING
HEALTHCARE
EQUIPMENT
PARTNER

Ultralight Carbon Fibre Wheeled Walker

Weighing just 5.7 kilograms, the stylish Lightweight Wheeled Walker is constructed from ultralight carbon fibre and designed for easy, every day use.

Height adjustable to best match the user, handle height is easily adjustable with an extension range of 120mm to suit users of average height and weight. Moulded handgrips with soft touch brakes are ideal for people with limited strength.

Features a soft, wide nylon seat and stable backrest for maximum comfort when seated. Also includes a large, detachable, weather resistant fabric bag that can hold up to 5kg of shopping or other personal belongings; and a stick holder to ensure safe storage and transportation of walking stick. This stylish European designer piece folds effortlessly with only one hand, simply pull the handle to collapse the walker for compact transportation or storage.

Features and Benefits

- Soft touch brake
- Stable backrest
- Soft wide seat
- Single handed folding mechanism
- 250mm front castors
- 200mm rear castors
- Concealed brake cables
- Reflectors for increased visibility in low light
- Detachable weather resistant fabric bag with reflective strip



Stands upright when folded



Easy fold for compact storage

Specifications

Height	865mm - 985mm
Length	690mm
Width	630mm
Seat Width	469mm
Seat Depth	210mm
Seat Height	540mm
Wheels	Front: 250mm Rear: 200mm
Folded Width	230mm
Weight	5.7kg
Bag & Backrest Weight	Additional 0.6kg
Bag SWL	5kg
SWL	136kg
Compliance	ARTG 176830 ISO 11199-2:2021
Warranty	12 months

rehab**hire**&sales

1 300 000 030

320 Lorimer Street Port Melbourne
contact@rehabhire.com.au
rehabhire.com.au

rehab**hire**&sales **REPAIRS AND SERVICE**

REGULAR MAINTENANCE OR REPAIRS

1 300 000 030

Unit 6 | 52 Wirraway Drive Port Melbourne
service@rehabhire.com.au
rehabhire.com.au/repairs-and-service

rehab**installation**

CAN MODIFY EVERY ROOM IN YOUR HOME TO
ENSURE IT IS ACCESSIBLE, SAFE, AND COMFORTABLE

1 300 000 030

320 Lorimer Street Port Melbourne
contact@rehabinstallation.com.au
rehabinstallation.com.au



Registered NDIS Provider