

rehabhire & sales

YOUR
LEADING
HEALTHCARE
EQUIPMENT
PARTNER

BetterLiving
by novis®

The Compact

Easy to lift, fold, and go



1300 000 030 | rehabhire.com.au

YOUR
LEADING
HEALTHCARE
EQUIPMENT
PARTNER

Lightweight Compact Wheeled Walker

Weighing just 7.2 kilograms, the stylish Compact Wheeled Walker is constructed from lightweight aluminium and designed for easy, every day use.

Height adjustable to best match the user, handle height is easily adjustable with an extension range of 120 mm to suit users of average height and weight. Moulded handgrips with soft touch brakes are ideal for people with limited strength.

Features a soft, wide nylon seat and stable backrest for maximum comfort when seated. Also includes a large, detachable, weatherproof fabric bag (with removable shoulder strap) that can hold up to 5 kg of shopping or other personal belongings; and a stick holder to ensure safe storage and transportation of walking stick. This stylish designer piece folds effortlessly with only one hand – simply pull the handle to collapse the walker, secure the seat clasp and then fold the handles down towards the wheels for double folding. The wheeled walker easily folds and can be carried without any extra securing needed, making it ideal for transport and storage.

Features and Benefits

- Soft touch brake
- Stable backrest
- Soft wide seat
- Single handed folding mechanism
- 200mm castors
- Concealed brake cables
- Reflectors for increased visibility in low light
- Detachable weather resistant fabric bag



Stands upright when folded

Easy fold for compact storage

Specifications

Height	820mm - 940mm
Length	696mm
Width	620mm
Seat Width	458mm
Seat Depth	250mm
Seat Height	545mm
Wheels	200mm
Folded Width	230mm
Weight	7.2kg
Bag & Backrest Weight	Additional 0.4kg
Bag SWL	5kg
SWL	136kg
Compliance	ARTG 176830 ISO 11199-2:2021
Warranty	12 months

rehabhire & sales

1 300 000 030

320 Lorimer Street Port Melbourne
contact@rehabhire.com.au
rehabhire.com.au

rehabhire & sales **REPAIRS AND SERVICE**

REGULAR MAINTENANCE OR REPAIRS

1 300 000 030

Unit 6 | 52 Wirraway Drive Port Melbourne
service@rehabhire.com.au
rehabhire.com.au/repairs-and-service

rehabinstallation

CAN MODIFY EVERY ROOM IN YOUR HOME TO
ENSURE IT IS ACCESSIBLE, SAFE, AND COMFORTABLE

1 300 000 030

320 Lorimer Street Port Melbourne
contact@rehabinstallation.com.au
rehabinstallation.com.au



Registered NDIS Provider