

rehabhire & sales

YOUR
LEADING
HEALTHCARE
EQUIPMENT
PARTNER

BetterLiving
by novis®

The Euro

Sleek design meets everyday functionality



1300 000 030 | rehabhire.com.au

YOUR
LEADING
HEALTHCARE
EQUIPMENT
PARTNER

Euro Lightweight Wheeled Walker

Weighing just 6.8 kilograms, the stylish Euro Lightweight Wheeled Walker is constructed from lightweight aluminium and designed for easy, every day use.

Height adjustable to best match the user, handle height is easily adjustable with an extension range of 100mm to suit users of average height and weight. Moulded handgrips with soft touch brakes are ideal for people with limited strength.

Features a soft, wide nylon seat and stable backrest for maximum comfort when seated. Also includes a large, detachable, weather resistant fabric bag (with removable shoulder strap) that can hold up to 4.5kg of shopping or other personal belongings; and a stick holder to ensure safe storage and transportation of walking stick. This stylish European designer piece folds effortlessly with only one hand, simply pull the handle to collapse the walker for compact transportation or storage.

Features and Benefits

- Soft touch brake
- Stable backrest
- Soft wide seat
- Single handed folding mechanism
- 200mm rear castors
- Concealed brake cables
- Reflectors for increased visibility in low light
- Detachable weather resistant fabric bag



Easy fold for compact storage



Stands upright when folded

Specifications

Height	815mm - 945mm
Length	780mm
Width	615mm
Seat Width	450mm
Seat Depth	250mm
Seat Height	540mm
Wheels	200mm
Folded Width	270mm
Weight	6.8kg
Bag SWL	4.5kg
SWL	136kg
SWL	136kg
Compliance	ARTG 176830 ISO 11199-2:2021
Warranty	12 months

rehab**hire**&sales

1300 000 030

320 Lorimer Street Port Melbourne
contact@rehabhire.com.au
rehabhire.com.au

rehab**hire**&sales **REPAIRS AND SERVICE**

REGULAR MAINTENANCE OR REPAIRS

1300 000 030

Unit 6 | 52 Wirraway Drive Port Melbourne
service@rehabhire.com.au
rehabhire.com.au/repairs-and-service

rehab**installation**

CAN MODIFY EVERY ROOM IN YOUR HOME TO
ENSURE IT IS ACCESSIBLE, SAFE, AND COMFORTABLE

1300 000 030

320 Lorimer Street Port Melbourne
contact@rehabinstallation.com.au
rehabinstallation.com.au



Registered NDIS Provider