

rehabhire & sales

YOUR  
LEADING  
HEALTHCARE  
EQUIPMENT  
PARTNER

BetterLiving  
by novis®

# The Premium

Comfort and class in every step



1300 000 030 | [rehabhire.com.au](http://rehabhire.com.au)

YOUR  
LEADING  
HEALTHCARE  
EQUIPMENT  
PARTNER

## Premium Lightweight Wheeled Walker

Weighing just 6.5 kilograms, the stylish Lightweight Wheeled Walker is constructed from lightweight aluminium and designed for easy, every day use.

Height adjustable to best match the user, handle height is easily adjustable with an extension range of 120mm to suit users of average height and weight. Moulded handgrips with soft touch brakes are ideal for people with limited strength.

Features a soft, wide nylon seat and stable backrest for maximum comfort when seated. Also includes a large, detachable, weather resistant fabric bag (with removable shoulder strap) that can hold up to 5kg of shopping or other personal belongings; and a stick holder to ensure safe storage and transportation of walking stick. This stylish designer piece folds effortlessly with only one hand, simply pull the handle to collapse the walker for compact transportation or storage.

### Features and Benefits

- Soft touch brake
- Stable backrest
- Soft wide seat
- Single handed folding mechanism
- 250mm front castors
- 200mm rear castors
- Concealed brake cables
- Reflectors for increased visibility in low light
- Detachable weather resistant fabric bag



Easy fold for compact storage

Stands upright when folded

### Specifications

Height	880mm - 970mm
Length	705mm
Width	610mm
Seat Width	461mm
Seat Depth	250mm
Seat Height	545mm
Wheels	Front: 250mm   Rear: 200mm
Folded Width	230mm
Weight	6.5kg
Bag & Backrest Weight	Additional 0.7kg
Bag SWL	5kg
SWL	136kg
Compliance	ARTG 176830   ISO 11199-2:2021
Warranty	12 months



# rehabhire & sales

1 300 000 030

320 Lorimer Street Port Melbourne  
contact@rehabhire.com.au  
rehabhire.com.au

# rehabhire & sales **REPAIRS AND SERVICE**

REGULAR MAINTENANCE OR REPAIRS

1 300 000 030

Unit 6 | 52 Wirraway Drive Port Melbourne  
service@rehabhire.com.au  
rehabhire.com.au/repairs-and-service

# rehabinstallation

CAN MODIFY EVERY ROOM IN YOUR HOME TO  
ENSURE IT IS ACCESSIBLE, SAFE, AND COMFORTABLE

1 300 000 030

320 Lorimer Street Port Melbourne  
contact@rehabinstallation.com.au  
rehabinstallation.com.au



Registered NDIS Provider