

rehabhire & sales

YOUR
LEADING
HEALTHCARE
EQUIPMENT
PARTNER



XXL Chair



1300 000 030 | rehabhire.com.au

YOUR
LEADING
HEALTHCARE
EQUIPMENT
PARTNER

The Cobi Rehab XXL chair is a strongly engineered and functional seating solution for bariatric users, designed to support a very high user weight of 325kg. The user-friendly design makes it easier to get in and out of the chair, providing greater independence and confidence. For added comfort, the height and seat depth of this chair can be easily adjusted to suit each individual. Designed to be easy to move, the chair uses high-tech, lightweight tubing, providing strength at the lowest weight.

Features and Benefits

- High user weight
- Height adjustable seat
- Butterfly armrests angle outwards for easier entry and exit
- Seat depth adjustable +/-60mm
- Long armrests protrude forward beyond the front of the seat for support
- Light-weight tubing to provide the highest strength at the lowest weight
- Polyurethane foam padded seat, backrest and armrests
- Easy to clean with a removable seat pad



Specifications

SKU	9544	9545
Height	790mm - 940mm	790mm - 940mm
Width	690mm	790mm
Depth	800mm - 1010mm	800mm - 1010mm
Seat Width	610mm	710mm
Seat Depth	400mm - 520mm	400mm - 520mm
Seat Height	440mm - 590mm	440mm - 590mm
Backrest Height	370mm	370mm
Width Between Arms	635mm	735mm
Frame	Powder-lacquered chromium molybdenum steel	
Seat Back Armrest	Polyurethane	
Product Weight	16.6kg	17.6kg
SWL	325kg	325kg
Compliance	DS EN 12182:2000 ISO 09 33 03	
Warranty	12 months limited	

rehabhire & sales

1 300 000 030

320 Lorimer Street Port Melbourne
contact@rehabhire.com.au
rehabhire.com.au

rehabhire & sales **REPAIRS AND SERVICE**

REGULAR MAINTENANCE OR REPAIRS

1 300 000 030

Unit 6 | 52 Wirraway Drive Port Melbourne
service@rehabhire.com.au
rehabhire.com.au/repairs-and-service

rehabinstallation

CAN MODIFY EVERY ROOM IN YOUR HOME TO
ENSURE IT IS ACCESSIBLE, SAFE, AND COMFORTABLE

1 300 000 030

320 Lorimer Street Port Melbourne
contact@rehabinstallation.com.au
rehabinstallation.com.au



Registered NDIS Provider