

rehabhire & sales

YOUR
LEADING
HEALTHCARE
EQUIPMENT
PARTNER



XXL Clean



1300 000 030 | rehabhire.com.au

YOUR
LEADING
HEALTHCARE
EQUIPMENT
PARTNER

The Cobi Rehab XXL Clean is a multi-purpose commode chair for bedside and shower use. Specifically designed for bariatric people, it supports a high user weight of 325kg. The user friendly design makes it simple to get in and out of the commode, enabling the user to maintain independence. For added comfort, the height and seat depth can be adjusted. The Clean commode is easy to move, and made of high-tech lightweight tubing that provides the highest strength at the lowest weight. Features and Benefits

- Seat depth adjustable +/-60mm
- Polyurethane foam padded seat, backrest and armrests
- Easy to clean with a removable seat pad
- Multipurpose commode chair that can be used as a bedside commode or a shower chair
- Comes with commode pan
- Height adjustable seat
- Butterfly armrests for easier entry and exit
- Long armrests protude forward beyond the front of the seat for support
- Light-weight tubing to provide the highest strength at the lowest weight



Specifications

SKU	9542	9543
Height	790mm - 940mm	790mm - 940mm
Width	690mm	790mm
Depth	800mm - 1010mm	800mm - 1010mm
Seat Width	610mm	710mm
Seat Depth	400mm - 520mm	400mm - 520mm
Seat Height	440mm - 590mm	440mm - 590mm
Backrest Height	370mm	370mm
Armrest Height	250mm	250mm
Distance Between Arms	690mm	790mm
Toilet Opening	Width: 240mm Length: 350mm	
Frame	Powder-lacquered chromium molybdenum steel	
Seat Back Armrest	Polyurethane	
Product Weight	14kg	15kg
SWL	325kg	325kg
Compliance	ARTG 34936 ISO 09 33 03	
Warranty	12 months limited	



Split Toilet Seat



Seat with Lid

rehabhire & sales

1 300 000 030

320 Lorimer Street Port Melbourne
contact@rehabhire.com.au
rehabhire.com.au

rehabhire & sales **REPAIRS AND SERVICE**

REGULAR MAINTENANCE OR REPAIRS

1 300 000 030

Unit 6 | 52 Wirraway Drive Port Melbourne
service@rehabhire.com.au
rehabhire.com.au/repairs-and-service

rehabinstallation

CAN MODIFY EVERY ROOM IN YOUR HOME TO
ENSURE IT IS ACCESSIBLE, SAFE, AND COMFORTABLE

1 300 000 030

320 Lorimer Street Port Melbourne
contact@rehabinstallation.com.au
rehabinstallation.com.au



Registered NDIS Provider