

rehabhire & sales

YOUR
LEADING
HEALTHCARE
EQUIPMENT
PARTNER



XXL Rollator King



1300 000 030 | rehabhire.com.au

YOUR
LEADING
HEALTHCARE
EQUIPMENT
PARTNER

The Cobi Rehab XXL Rollator King is a reinforced rollator designed for bariatric users. With a lightweight and discrete design, the rollator supports a very high user weight of 325kg. The frame fits through doorways and folds easily for transport and storage by lifting the seat. Often a critical point of rollators for bariatric users is the effectiveness of the brakes. The 510mm wide seat has a back arch that provides comfort and safety, preventing the user's weight from shifting too far back.

Features and Benefits

- Iron forks and four large 200mm wheels ensure a smooth drive and safe handling
- Feather-touch hand brakes are highly responsive for superior control and support
- Folds down for transport and storage
- Height adjustable in 25mm intervals
- Powder-lacquered aluminium frame
- Large supportive seat with back support
- Convenient basket for personal belongings
- Reflectors fitted to handles for safety

Specifications

SKU	9555	9661
Height	775mm - 990mm	775mm - 990mm
Width	740mm	840mm
Length	850mm	850mm
Seat Width	510mm	510mm
Seat Depth	370mm	370mm
Seat Height	510mm	510mm
Handle Height	775mm - 990mm	775mm - 990mm
Width Between Handles	610mm	710mm
Wheels	200mm	
Frame	Steel SAE4130	
Seat Back	Polyurethane	
Folded Dimensions	740mm	740mm
Product Weight	15.7kg	15.7kg
Basket SWL	5kg	5kg
SWL	325kg	325kg
Compliance	ARTG 37952 ISO 12 06 06 – 11	
Warranty	12 months limited	



rehab**hire** & sales

1 300 000 030

320 Lorimer Street Port Melbourne
contact@rehabhire.com.au
rehabhire.com.au

rehab**hire** & sales **REPAIRS AND SERVICE**

REGULAR MAINTENANCE OR REPAIRS

1 300 000 030

Unit 6 | 52 Wirraway Drive Port Melbourne
service@rehabhire.com.au
rehabhire.com.au/repairs-and-service

rehab**installation**

CAN MODIFY EVERY ROOM IN YOUR HOME TO
ENSURE IT IS ACCESSIBLE, SAFE, AND COMFORTABLE

1 300 000 030

320 Lorimer Street Port Melbourne
contact@rehabinstallation.com.au
rehabinstallation.com.au



Registered NDIS Provider