

rehabhire & sales



YOUR  
LEADING  
HEALTHCARE  
EQUIPMENT  
PARTNER

# XXL Shower Commode



1300 000 030 | [rehabhire.com.au](http://rehabhire.com.au)

YOUR  
LEADING  
HEALTHCARE  
EQUIPMENT  
PARTNER

The Cobi Rehab® XXL shower commode has been specifically designed as a showering solution for bariatric people, supporting a very high user weight of 325kg. The functional, user-friendly design provides good support for a mobile user to get into and out of the commode independently. The armrests are positioned to protrude out 70mm from the seat base, making them easy to access. The practical design of the commode ensures transfers, positioning and cleaning is made easier for the carer.

## Features and Benefits

- The extra 70mm between the edge of seat to the armrests and the 50mm high PU seat allows additional space for the user
- The footrests can support 100kg each and are covered with a tough PU foam for both durability and comfort
- Swing away footrests can be removable and width adjustable up to 1100mm apart
- Wide internal seat accommodates the user, while the external width will fit through most modern doors
- Large space between back rest and seat for space to fit a sling around the user
- Open back toilet seat with a long open hole provides easy access for cleaning and hygiene
- Armrests extend 70mm beyond the front of the seat base for support when transferring in and out
- Easily cleaned PU foam on toilet seat, backrest, armrest, and footrest cover can be replaced if required



Split Toilet Seat



Seat with Lid

## Specifications

SKU	9540	9541
Height	1000mm	1000mm
Width	680mm	780mm
Depth	1170mm	1170mm
Seat Width	610mm	710mm
Seat Depth	560mm	560mm
Seat Height	550mm	550mm
Backrest Height	480mm	480mm
Toilet Opening	Width: 240mm   Length: 350mm	
Footrest Height	400mm - 480mm	400mm - 480mm
Footrest Width	580mm - 1100mm	580mm - 1100mm
Wheels	125mm	125mm
Frame	Powder-lacquered chromium molybdenum steel	
Seat   Back   Armrest	Polyurethane	
Product Weight	27kg	30kg
SWL	325kg	325kg
Compliance	ARTG 34936   ISO 09 12 03	
Warranty	12 months limited	





# rehab**hire**&sales

1300 000 030

320 Lorimer Street Port Melbourne  
contact@rehabhire.com.au  
rehabhire.com.au

# rehab**hire**&sales **REPAIRS AND SERVICE**

REGULAR MAINTENANCE OR REPAIRS

1300 000 030

Unit 6 | 52 Wirraway Drive Port Melbourne  
service@rehabhire.com.au  
rehabhire.com.au/repairs-and-service

# rehab**installation**

CAN MODIFY EVERY ROOM IN YOUR HOME TO  
ENSURE IT IS ACCESSIBLE, SAFE, AND COMFORTABLE

1300 000 030

320 Lorimer Street Port Melbourne  
contact@rehabinstallation.com.au  
rehabinstallation.com.au



Registered NDIS Provider