

rehabhire & sales



YOUR
LEADING
HEALTHCARE
EQUIPMENT
PARTNER

XXL Stand Tall Walking Frame



1300 000 030 | rehabhire.com.au

YOUR
LEADING
HEALTHCARE
EQUIPMENT
PARTNER

The Cobi Rehab® XXL Stand Tall walking frame is a multifunctional mobility aid that can be used for walking and sit to stand transfers. The Stand Tall gives greater independence as it allows users to engage in their own transfers. Specifically designed for bariatric users, this mobility aid has a user weight of 325kg. 32mm tubing ensures the frame is durable and long lasting as well as easy to clean. The frame has been cleverly designed to provide exceptional support and to minimise the risk of injury to both the carer and the user. Two built-in hand grips at different heights, allow the user to rise, to a standing position in small stages. The brakes are easy to access and can be locked for greater control and stability. The frame height is adjustable and the width and angle of the forearm padding can also be adapted to suit each individual.

Features and Benefits

- Easily adjustable frame height
- The width and angle of the forearm padding is adjustable
- Lockable brakes
- Powder-lacquered chromium molybdenum steel
- Two PU foam hand grips on each side with multiple grip combinations
- Forearm padding
- Durable, easily cleaned frame constructed from 32mm diameter tubing



Specifications

Height	1000mm - 1300mm
Width	760mm
Length	730mm
Width between underarm support centre to centre	610mm - 810mm
Internal Width	720mm
Stand Tall Turning Radius	2000mm
Wheels	Front: 100mm Rear: 125mm
Frame	Powder-lacquered chromium molybdenum steel
Underarm Support Grip	Polyurethane
Product Weight	20kg
SWL	325kg
Compliance	ARTG 31117
Warranty	12 months limited

rehab**hire**&sales

1300 000 030

320 Lorimer Street Port Melbourne
contact@rehabhire.com.au
rehabhire.com.au

rehab**hire**&sales **REPAIRS AND SERVICE**

REGULAR MAINTENANCE OR REPAIRS

1300 000 030

Unit 6 | 52 Wirraway Drive Port Melbourne
service@rehabhire.com.au
rehabhire.com.au/repairs-and-service

rehab**installation**

CAN MODIFY EVERY ROOM IN YOUR HOME TO
ENSURE IT IS ACCESSIBLE, SAFE, AND COMFORTABLE

1300 000 030

320 Lorimer Street Port Melbourne
contact@rehabinstallation.com.au
rehabinstallation.com.au



Registered NDIS Provider