

rehabhire & sales

YOUR
LEADING
HEALTHCARE
EQUIPMENT
PARTNER



XXL Walking Frame



1300 000 030 | rehabhire.com.au

YOUR
LEADING
HEALTHCARE
EQUIPMENT
PARTNER

The Cobi Rehab XXL walking frame has been designed to offer bariatric users a reliable walking aid that promotes mobility. The Cobi Rehab XXL walking frame's strong, functional and lightweight design has been specifically developed to support a high user weight of 325kg. The frame height can be adjusted without tools from 760mm to 1010mm, and folds easily for transport and storage. The frame is available in two widths: 740mm and 840mm.

Features and Benefits

- Adjustable padded handles made with PU foam
- Rubber tips on front and rear
- Very high user weight for bariatric users
- Foldable for transport and storage
- Lightweight frame
- Height adjustable without tools

Specifications

SKU	9553	9554
Height	760mm - 1010mm	760mm - 1010mm
Width	740mm	840mm
Depth	500mm	500mm
Inside Width	540mm	640mm
Handle Height	760mm - 1010mm	760mm - 1010mm
Frame	Powder-lacquered light steel	
Product Weight	3kg	3kg
SWL	325kg	325kg
Compliance	ARTG 37938 ISO 09 12 03	
Warranty	12 months limited	



rehabhire & sales

1 300 000 030

320 Lorimer Street Port Melbourne
contact@rehabhire.com.au
rehabhire.com.au

rehabhire & sales **REPAIRS AND SERVICE**

REGULAR MAINTENANCE OR REPAIRS

1 300 000 030

Unit 6 | 52 Wirraway Drive Port Melbourne
service@rehabhire.com.au
rehabhire.com.au/repairs-and-service

rehabinstallation

CAN MODIFY EVERY ROOM IN YOUR HOME TO
ENSURE IT IS ACCESSIBLE, SAFE, AND COMFORTABLE

1 300 000 030

320 Lorimer Street Port Melbourne
contact@rehabinstallation.com.au
rehabinstallation.com.au



Registered NDIS Provider