



Booster

by Kalagon®

Booster is designed to reduce anxiety when using your existing air floatation cushion by maintaining your clinician's recommendation for pressure, and eliminating the need to manually adjust it. This device improves your existing seating without making big changes, and can be installed onto most full-sized single-cell air flotation wheelchair cushions.

Features and Benefits

- Automatic leak detection
- Auto-inflation upon leak detection
- Extended battery life
- Automatic adjustments for environmental changes and when you shift
- Personalised settings in the Kalagon® app based on your preferences and your clinician's guidance
- Checks status more than once a second



Specifications

Total Volume	Length: 118.036 mm Width: 57.32 mm Height (without clip):37.31 mm
Weight	Less than 0.2kg for the bare device without the clip
Hose Length	Delivered at a length 300mm which can be cut to preferred length
USB Charging	Chargers adhering to the USB-IF standard. Compatible with USB C charging plugs.
Battery Chemistry	Lithium-Polymer (LIP)
Electromagnetic Compatibility (EMC)	Classified as a Class 1, Rule 13 medical device per the EU MDR. EMC testing performed to IEC 60601-1-2 testing standards.
Air Travel and Battery Safety	UN38.3 AMS
SWL	Defined by third-party wheelchair cushion manufacturer

rehabhire & sales

1 300 000 030

320 Lorimer Street Port Melbourne
contact@rehabhire.com.au
rehabhire.com.au

rehabhire & sales **REPAIRS AND SERVICE**

REGULAR MAINTENANCE OR REPAIRS

1 300 000 030

Unit 6 | 52 Wirraway Drive Port Melbourne
service@rehabhire.com.au
rehabhire.com.au/repairs-and-service

rehabinstallation

CAN MODIFY EVERY ROOM IN YOUR HOME TO
ENSURE IT IS ACCESSIBLE, SAFE, AND COMFORTABLE

1 300 000 030

320 Lorimer Street Port Melbourne
contact@rehabinstallation.com.au
rehabinstallation.com.au



Registered NDIS Provider