

rehabhire&sales

YOUR
LEADING
HEALTHCARE
EQUIPMENT
PARTNER



HealthFlex Pivot



1300 000 030 | rehabhire.com.au

YOUR
LEADING
HEALTHCARE
EQUIPMENT
PARTNER

The innovative Oscar ‘Pivot’ mattress hinge is designed to significantly reduce skin ‘shear strain’ when raising and lowering your electric bed backrest.

By allowing the mattress surface to bend and glide freely over adjustable bed bases, a more comfortable user experience is achieved.

Reduced body movement across the mattress surface also reduces ‘sliding friction’ discomfort between the skin and mattress surface.

Features and Benefits

- Low profile 22cm mattress height
- Foam box edge support provides excellent transfer stability
- A clinical mattress designed to suit the aesthetics of your home environment
- ARTG number (Australian Registry of Therapeutic Goods) embroidered on the foot end border panel
- Mattress Law tag sewn at the head end of the mattress with clear instructions and user guidelines
- SWL tag

Joyce HygroFlex® Foam

The Australian Made, Oscar Health Flex PIVOT is available in three feels, soft, medium and firm. All feels include Joyce HygroFlex® foam tested by RMIT University to show the following benefits:

- More effective pressure redistribution than traditional memory foams
- Superior Airflow which helps diffuse body heat more rapidly
- Less ‘sinking feeling’ than experienced with most memory foams.
- Reduced surface shear and friction



- Sanitised stretch knit fabric.
- A super soft foam gives the quilting depth and a soft, luxurious feel.
- Premium high density Australian foams including HygroFlex open cell foam, which enhance comfort and provide excellent support.
- Solid high density foam foundation layer is built with firm edges, designed to suit articulating beds, and to help prevent users ‘bottoming out’ onto the bed base.



Specifications

Size	Total Length	Total Width	Total Height	SWL
Single	1880mm	920mm	220mm	150kg
Long Single	2030mm	920mm	220mm	150kg
King Single	2030mm	1060mm	220mm	150kg
Double	1870mm	1370mm	220mm	250kg
Long Double	2030mm	1370mm	220mm	250kg
Queen	2030mm	1530mm	220mm	250kg
Split Queen	2030mm	760mm	220mm	150kg per side
Compliance	ARTG 448347			
Warranty	Foams: 10 years Cover: 5 years			



Notes

Lined area for notes, consisting of 25 horizontal lines.

rehab**hire**&sales

1 300 000 030

320 Lorimer Street Port Melbourne
contact@rehabhire.com.au
rehabhire.com.au

rehab**hire**&sales **REPAIRS AND SERVICE**

REGULAR MAINTENANCE OR REPAIRS

1 300 000 030

Unit 6 | 52 Wirraway Drive Port Melbourne
service@rehabhire.com.au
rehabhire.com.au/repairs-and-service

rehab**installation**

CAN MODIFY EVERY ROOM IN YOUR HOME TO
ENSURE IT IS ACCESSIBLE, SAFE, AND COMFORTABLE

1 300 000 030

320 Lorimer Street Port Melbourne
contact@rehabinstallation.com.au
rehabinstallation.com.au



Registered NDIS Provider