

rehabhire & sales

YOUR  
LEADING  
HEALTHCARE  
EQUIPMENT  
PARTNER



QUANTUM®  
**EDGE 4**



1300 000 030 | [rehabhire.com.au](http://rehabhire.com.au)

YOUR  
LEADING  
HEALTHCARE  
EQUIPMENT  
PARTNER



## The Most Advanced Edge Yet

The Edge® 4 is our most advanced model yet. With 178mm front and rear casters and powerful 4-pole motors, it delivers speeds up to 10km/h and confidently tackles everyday obstacles. A new battery tray allows easier access and supports 75 Ah batteries for extended range. LED accent lights are integrated into the casters and fenders for greater visibility and style. Enhanced suspension and damping smooth out rough terrain, giving you comfort, control, and performance wherever you go.















### Features and Benefits

- Independent SRS (Smooth Ride Suspension)
- Standard led front castor & rear fender lighting package
- USB charger included
- Wide array of color choices
- Q-Logic 3 Drive Controls
- Made in Australia of foreign and domestic parts

### Options

- TRU-Balance® 4 Power positioning:
- 10° or 30° anterior tilt Options with memory seating
- TRU-Balance® 3 Power Positioning:
- iLevel® option allows 30.5 cm of seat elevation while at walking speed (7.6km/h) Power Tilt, Recline, and Elevating Legrests that can power to floor for easy transfers.
- Accu-Trac advanced tracking technology
- NEW backup camera
- Additional options available on scripts / configurator

### Colours

	Forest Green		Back in Black
	Cosmic Purple		Night Train
	Bermuda Blue		Steel Blue
	Lemon Crush		Royal Blue
	Evolution Orange		Inferno Red
	Blue Jean		Ice Pink
	Pretty in Pink		Ghost Gray



Front Castor Accent Lights



Rear Fender Accent Lights



Upgraded Suspension



178mm Castors



New Battery Tray



USB Charger

## Specifications

Drive Wheels	356mm pneumatic or flat free tyres Optional: 356mm Community flat free tyres available
Castor Wheels	178mm flat free
Maximum Speed <sup>1</sup>	<10km/h <7.25km/h with iLevel
Range per charge <sup>1,8</sup>	29.2 km with Group 24 batteries
Ground Clearance	73mm
Turning Radius	579mm
Length <sup>2</sup>	972.06mm
Base Width	610.11mm
<b>Seating Sizes</b>	
TRU Balance® 4 Anterior Tilt <sup>5</sup>	Width: 400mm - 558mm Depth: 355mm - 558mm
TRU Balance® 3 Static Seat	Width: 305mm - 558mm Depth: 305mm - 508mm
TRU Balance® 3 Power Tilt	Width: 305mm - 558mm Depth: 305mm - 558mm
TRU Balance® 3 Tilt and iLevel®	Width: 305mm - 558mm Depth: 305mm - 558mm
<b>Seat to Floor Heights</b>	
TRU Balance® 4 Anterior Tilt <sup>5</sup>	463mm - 539mm
TRU Balance® 3 Power Tilt <sup>5</sup>	450mm - 501mm
TRU Balance® 3 Tilt and iLevel® <sup>5</sup>	444.5mm - 495mm
<b>Battery Size <sup>3</sup></b>	
Battery Weight <sup>7</sup>	MK 70Ah (2) 23.3kg each
Available Electronics <sup>6</sup>	120A Q-Logic 3
Battery Charger	8A Off-board
Motor Packages	4-Pole
Base Weight exc. Batteries	70kg
Weight Capacity <sup>4</sup>	150kg

<sup>1</sup> Speed varies with user weight, terrain type, battery charge, battery condition and tyre pressure.

<sup>2</sup> Without front riggings.

<sup>3</sup> Battery capacity may change when adding power positioning systems.

<sup>4</sup> Weight capacity is dictated by the selection of seating.

<sup>5</sup> All measurements +/- 3.18 with 457.2mm x 457.2mm seat, fixed towers, 150kg. weight capacity and solid tyres (add 9.52 with pneumatic tyres). Widths may increase seat-to-floor height. Contact Customer Service for details.  
Listed measurements are in the lowest position; does not include elevating system unless specified.

<sup>6</sup> All Quantum Rehab® units can be ordered with specialty controls.

<sup>7</sup> Battery weight may vary +/- 1.5kg based on manufacturer.

<sup>8</sup> Range shown is measured using Group 24 batteries.

<sup>9</sup> Pneumatic not available in community tyre option.



# rehab**hire**&sales

1 300 000 030

320 Lorimer Street Port Melbourne  
contact@rehabhire.com.au  
rehabhire.com.au

# rehab**hire**&sales **REPAIRS AND SERVICE**

REGULAR MAINTENANCE OR REPAIRS

1 300 000 030

Unit 6 | 52 Wirraway Drive Port Melbourne  
service@rehabhire.com.au  
rehabhire.com.au/repairs-and-service

# rehab**installation**

CAN MODIFY EVERY ROOM IN YOUR HOME TO  
ENSURE IT IS ACCESSIBLE, SAFE, AND COMFORTABLE

1 300 000 030

320 Lorimer Street Port Melbourne  
contact@rehabinstallation.com.au  
rehabinstallation.com.au



Registered NDIS Provider