

rehabhire & sales

YOUR
LEADING
HEALTHCARE
EQUIPMENT
PARTNER



M5

Bariatric Lift Recliner



1300 000 030 | rehabhire.com.au

YOUR
LEADING
HEALTHCARE
EQUIPMENT
PARTNER

The M5 is the first choice for a premium bariatric lift chair, with its class leading Australian made four motor mechanism, culmination of years of research, testing of design and development by our in house team of engineers. Reliability was a key focus of the M5 development process which is why we put it through thousands of lift cycles while weighted above its rated capacity at each stage of development. This has resulted in the most durable bariatric mechanism on the market.

Features and Benefits

- 4 medical-grade electric motors
- Separate control of the backrest, lift, and footrest
- Right side controller
- 24V rechargeable battery backup
- Memory foam seat cushion
- Pocket spring suspension
- 75mm heavy-duty castors
- Armrest covers
- Showwood arms
- Available in a wide range of Versotela fabrics

Modification Options

- Cushion cutout
- Memory foam back, armrests, and footrest
- Left side controller
- Gluteal recess
- No showwood



All VersoTela fabrics come standard with a suite of high-end specifications, such as, antimicrobial coating, mold resistance, scratch resistance, and fire resistance. The wide range of collections each boast unique properties that are carefully designed to suit every customer's needs. Whether it be built for use in a care facility, or as a functional, but stylish addition to the living room, the VersoTela fabric options are sure to tick all the boxes when it comes to crafting your perfect chair.



Specifications

| Size | 650 | 650/70 | 850 | 850/70 |
|---------------------------|-------|--------|-------|--------|
| Seat Height | 500mm | 570mm | 500mm | 570mm |
| Seat Width | 650mm | 650mm | 850mm | 850mm |
| Seat Depth | 530mm | 600mm | 530mm | 600mm |
| Seat Depth Min. Adj. | -50mm | -50mm | -50mm | -50mm |
| Seat Depth Max. Adj. | +50mm | +50mm | +50mm | +50mm |
| Backrest Height | 636mm | 647mm | 688mm | 718mm |
| Backrest Height Min. Adj. | -25mm | -25mm | -25mm | -25mm |
| Backrest Height Max. Adj. | +75mm | +75mm | +75mm | +75mm |
| Arm Height | 190mm | 190mm | 190mm | 190mm |
| Arm Height Min. Adj. | -25mm | -25mm | -25mm | -25mm |
| Footrest Length | 950mm | 1020mm | 950mm | 1020mm |
| Footrest SWL | 80kg | 80kg | 80kg | 80kg |
| SWL | 320kg | 320kg | 320kg | 320kg |

Warranty

Proof of Purchase:

This Oscar Chair is fitted with a compliance card on the rear bottom rail. All cards have a unique ID number called the "SO" number. Warranty is voided if this card is removed.

Electronics and Handset:

2 year warranty on all electrical.

Foam/Frame:

5 year warranty on material and workmanship, excluding general wear and tear.

Fabric:

2 year warranty on the Versotela fabric ranges when the care labels have been followed

Battery Pack:

2 year warranty on the battery pack.

Terms

All warranties apply from the date of sale. Where practical the Rehab Hire & Sales will fit warranty replacement parts to minimise repair times for the end user. All warranty claims are to be raised with Rehab Hire & Sales. Rehab Hire & Sales will liaise with the manufacturer to process the claim.

Exclusions

- Damage due to incorrect use of the chair.
- General wear and tear.
- Weight capacity having been exceeded.
- Any modification made to the product by the consumer or a 3rd party without prior written approval from Rehab Hire & Sales to the nature and extent of the modification.
- Minor variations in the finish of the product which does not affect the products value of fitness for its intended purpose.

rehabhire & sales

1 300 000 030

320 Lorimer Street Port Melbourne
contact@rehabhire.com.au
rehabhire.com.au

rehabhire & sales **REPAIRS AND SERVICE**

REGULAR MAINTENANCE OR REPAIRS

1 300 000 030

Unit 6 | 52 Wirraway Drive Port Melbourne
service@rehabhire.com.au
rehabhire.com.au/repairs-and-service

rehabinstallation

CAN MODIFY EVERY ROOM IN YOUR HOME TO
ENSURE IT IS ACCESSIBLE, SAFE, AND COMFORTABLE

1 300 000 030

320 Lorimer Street Port Melbourne
contact@rehabinstallation.com.au
rehabinstallation.com.au



Registered NDIS Provider