

rehabhire & sales

YOUR  
LEADING  
HEALTHCARE  
EQUIPMENT  
PARTNER

ottobock.

# Pointer Sports Chair



1300 000 030 | [rehabhire.com.au](http://rehabhire.com.au)

YOUR  
LEADING  
HEALTHCARE  
EQUIPMENT  
PARTNER

## The Pointer Sports Chair Range.

The Pointer sports chair offers a comprehensive range of settings and options so that your needs as a sportsperson can be fully met. The seat height, lower leg length and centre of gravity can all be individually adjusted to meet your requirements and body measurements.

### Benefits at a glance

- Configurable sports wheelchair
- Seat angle, centre of gravity and back angle with different setting options
- Adaptable back upholstery



Adaptable, height and angle-adjustable leg support



Adjustable centre of gravity



Adjustable back support angle



Individually adjustable seat height



### Specifications

Seat Width	340mm - 460mm
Seat Depth	360mm - 480mm
Seat Height   Anterior	Depth: 450mm   470mm   490mm Height: 490mm   510mm   530mm
Seat Height   Posterior	Depth: 370mm - 450mm Height: 430mm - 530mm
Back Height	250mm - 500mm
Lower Leg Length	340mm - 470mm
Front Frame Offset	20mm
Wheel Camber	18°
SWL	120kg



# rehabhire & sales

1 300 000 030

320 Lorimer Street Port Melbourne  
contact@rehabhire.com.au  
rehabhire.com.au

# rehabhire & sales **REPAIRS AND SERVICE**

REGULAR MAINTENANCE OR REPAIRS

1 300 000 030

Unit 6 | 52 Wirraway Drive Port Melbourne  
service@rehabhire.com.au  
rehabhire.com.au/repairs-and-service

# rehabinstallation

CAN MODIFY EVERY ROOM IN YOUR HOME TO  
ENSURE IT IS ACCESSIBLE, SAFE, AND COMFORTABLE

1 300 000 030

320 Lorimer Street Port Melbourne  
contact@rehabinstallation.com.au  
rehabinstallation.com.au



Registered NDIS Provider