

rehabhire & sales

showerbuddy®

YOUR  
LEADING  
HEALTHCARE  
EQUIPMENT  
PARTNER

SB1



1300 000 030 | [rehabhire.com.au](http://rehabhire.com.au)

YOUR  
LEADING  
HEALTHCARE  
EQUIPMENT  
PARTNER

## Transfer to a shower with a tray step

The ShowerBuddy SB1 shower transfer system will turn your existing bathroom into an accessible bathroom allowing you to enjoy showering again without the excessive cost, delays and stress of a full bathroom remodel. If your shower has a step or upstand preventing roll-in access, our ShowerBuddy SB1 will easily slide you, in a seated position, across the step and into the shower. Once transferred into the shower the bridge and wheeled base can be disconnected and rolled away, allowing the shower door to close for a more enjoyable shower experience. This chair can be rotated when inside the shower allowing you to face the water outlet or for an assistant to help with cleaning.

## Features and Benefits

- Chair and bridge design
- Rolls over most standard toilets
- Low maintenance
- Removable, lockable, flipback arms
- Height adjustable
- Safety belts: lap and chest belts
- Locking castor wheels: 125mm x 4
- Safety mechanisms
- No rust aluminium frame | Grade 6061 T6
- Stainless steel parts and hardware | Grade 304
- Rotates 360° inside the shower (subject to space)
- Fold up footrest
- Commode system
- Tool-free assembly
- Set up and use in under an hour
- No heavy lifting at shower time
- No bathroom remodelling required





# rehabhire & sales

1 300 000 030

320 Lorimer Street Port Melbourne  
contact@rehabhire.com.au  
rehabhire.com.au

# rehabhire & sales **REPAIRS AND SERVICE**

REGULAR MAINTENANCE OR REPAIRS

1 300 000 030

Unit 6 | 52 Wirraway Drive Port Melbourne  
service@rehabhire.com.au  
rehabhire.com.au/repairs-and-service

# rehabinstallation

CAN MODIFY EVERY ROOM IN YOUR HOME TO  
ENSURE IT IS ACCESSIBLE, SAFE, AND COMFORTABLE

1 300 000 030

320 Lorimer Street Port Melbourne  
contact@rehabinstallation.com.au  
rehabinstallation.com.au



Registered NDIS Provider