

rehabhire & sales

PRIDE[®]
MOBILITY

YOUR
LEADING
HEALTHCARE
EQUIPMENT
PARTNER

PATHRIDER[®] **10**
DELUXE

EZ TURN
TECHNOLOGY



1300 000 030 | rehabhire.com.au

Style, Comfort and Performance

The Pathrider® 10 Deluxe scooter offers a sleek and sporty design with smooth, reliable performance. Equipped with advanced standard features like EZ turn technology, rugged suspension, a front and rear lighting system for high visibility, and an easy-to-use wraparound tiller, the Pathrider® 10 Deluxe delivers exceptional value and convenience.

Specifications

Width ³	562mm
Length ³	1194mm
Seat Width	442mm
Seat Depth	430mm
Seat to Deck Height	438mm – 476mm
Seat to Ground Height	600mm – 638mm
Ground Clearance ⁵	825mm at frame
Turning Radius	1155mm
Tyres Front	264.16mm
Tyres Rear	273.05mm
Range per Charge ¹	<25km
Maximum Speed ¹	8.9kph
Battery Charger	Off-board, 3.5A (32Ah) or Off-board, 5A (40Ah)
Battery Requirements ⁴	(2) 31Ah (standard) or (2) 40Ah (optional extra)
Battery Weight	10.5kg (32Ah) or 14kg (40Ah)
Weight without Battery ²	84.8kg
Colour	Viper Blue
SWL	158.7kg
Compliance	AS3695.2:2013
Warranty	2 years Batteries: 12 months

1. Varies with user weight, terrain type, battery amp hour (Ah), battery charge, battery condition and tyre condition.
2. Includes standard seat weight.
3. Width of base (not including seat).
4. Deep cycle (AGM or gel cell type recommended)
5. Does not account for anti-tip clearance.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tyre pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of 3%.



Front and Rear
Suspension



High Back Seat with
Swivel and Slide



Full LED Lighting
Package



1300 000 030

320 Lorimer Street Port Melbourne
contact@rehabhire.com.au
rehabhire.com.au



1300 000 030

Unit 6 | 52 Wirraway Drive Port Melbourne
service@rehabhire.com.au
rehabhire.com.au/repairs-and-service



1300 000 030

320 Lorimer Street Port Melbourne
contact@rehabinstallation.com.au
rehabinstallation.com.au