

YOUR
LEADING
HEALTHCARE
EQUIPMENT
PARTNER

Extremely user-friendly and high performance skin protection

- High performance skin protection by even pressure re-distribution
- Extreme ease of use
- Supporting the pelvis and thighs for a stable sitting posture



Colour coded label
Easily shows if cushion and cover are the same size

Safe to use in all directions: top = bottom; left = right; front = back.



PRESSURE INJURY TREATMENT PLAN

SOLID SHAPE

READY FOR USE

100% BREATHABLE

MACHINE WASHABLE

Product performance

Skin protection / Anti decubitus	
Pressure redistribution	◀◀◀◀◀
Shear force reduction	◀◀◀◀◀◀
Microclimate control	◀◀◀◀◀
Impact & vibration damping	◀◀◀◀◀

Positioning	
Lateral stability	◀◀◀◀◀◀
Frontal stability	◀◀◀◀◀
Sliding reduction	◀◀◀◀◀
Allows manual adjustments	N/A
Increasing endurance	◀◀◀◀◀◀

Maintenance & Safety	
Durability	◀◀◀◀◀
Reliability	◀◀◀◀◀
Functionality when misused	◀◀◀◀◀
Easy to maintain	◀◀◀◀◀
Level of hygiene	◀◀◀◀◀

Functionality & application	
Comfort	◀◀◀◀◀
Ultra lightweight	◀◀◀◀◀
Allows easy transfer	◀◀◀◀◀◀
Available standard sized	◀◀◀◀◀◀
Maximum user weight	◀◀◀◀◀

Technical specifications

Thickness	8 cm
Average weight	800 grams
Max. user weight	250 kg
Machine washable	60°C max.

Cushion sizes

- Standard
- Not available

Inch	width				
	16	17	18	19	20
CM	40	43	45	48	50
16	40				
17	43	43			
18	45		45		
19	48			48	
20	50				50

Available covers

Comfair Cover and Incotec Cover.

The Vicair 4 consists of 4 identical individual compartments designed to contour under loading. The cushion is always placed correctly, as all four sides can be either front or back even the top and bottom are identical. The configuration of the four compartments is designed to work together to strike a unique balance between comfort, preserving skin and underlying tissue and pelvic support. The four identical, unadjustable, individual compartments allow for automatic positioning on the seating surface. Vicair 4 cushions and covers have the unique feature of colour coded labels, easily showing if the cushion and cover are the same size. This is considered to be beneficial in environments where many wheelchair cushions are in use at once, such as care centres.



1300 000 030

320 Lorimer Street Port Melbourne
contact@rehabhire.com.au
rehabhire.com.au



1300 000 030

Unit 6 | 52 Wirraway Drive Port Melbourne
service@rehabhire.com.au
rehabhire.com.au/repairs-and-service



1300 000 030

320 Lorimer Street Port Melbourne
contact@rehabinstallation.com.au
rehabinstallation.com.au